Serie **SILVER+** TG-1200

KARTSANA

Life-moving innovation

C/ Narcís Monturiol, 34
08192 Sant Quirze del Vallès
Barcelona (España)
+34 93 715 86 72
info@kartsana.com
www.kartsana.com

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Electric stretcher

Model TG-1200

Please read these maintenance instructions carefully before using the product and keep them for future reference.

Electric

CONTENTS

1. INTRODUCTION	4	5. SPARE PARTS	28
1.1 Foreword	4	6. ACCESSORIES	32
1.2 Intended product use	5	7. SYSTEM ERRORS	33
1.3 Liability and warranty	5	7.1 Error codes	33
1.4 Specifications	5	7.2 Troubleshooting	33
1.5 Attention	5	7.3 Motor protection	35
1.6 Summary of safety precautions	5	7.4 Ambulance entry	35
2. INSTRUCTION MANUAL	10	7.5 Battery charging	36
2.1 Technical characteristics of the stretcher	10		
2.2 Vehicle configuration	12		
2.3 Operation and handling of the stretcher	13		
2.3.1 Fitting the batteries	13		
2.3.2 Setting up and control of the stretcher	14		
2.3.3 Adjustable head	15		
2.3.4 Adjustable leg	15		
2.3.5 Safety railing	16		
2.3.6 Drip holder	17		
2.3.7 Safety belts	18		
2.3.8 Brake lever for rear wheels	19		
2.3.9 Lever for locking the front wheels			
2.3.10 Mounting of the stretcher on the R-1200 rail	19		
2.3.11 Transferring the patient to the stretcher	20		
	23		
3. FAULTS	24		
4. CARE AND MAINTENANCE	25		
4.1 Contact information	25		
4.2 Environmental conditions	25		
4.3 Cleaning	25		
4.4 Maintenance	26		
4.4.1 Precautionary maintenance	26		
4.4.2 Service maintenance	27		
4.5 Transportation and storage	27		
4.6 Disposal	27		

1. INTRODUCTION

1.1 Foreword

The **SILVER + (TG-1200)** stretcher has been specifically designed for patient transport. The product complies with the MDR (EU) Directive 2017/745. It is an electro-hydraulic model that enables the adjustment of the working height according to individual need, with just two push buttons. There is the option of manual operation for possible malfunctions of the automatic system.

The SILVER stretcher model must be attached to the ambulance by means of the **SILVER+ R-1200** rail.

This stretcher has been tested in accordance with the UNE-EN 1789 standard.

All information on treatment, disinfection and maintenance is based on our current experience and understanding.

We reserve the right to make variations in the technical specifications of the stretcher in order to improve the product.

SILVER+ information sticker

Life-moving innovation

KARTSANA, S.L.
C/Narcís Monturiol 34
08192 St. Quirze del Vallès
Barcelona (Spain)

Model: **TG-1200**
Commercial name: **SILVER+**
Voltage: **24V DC**
Max. Power Input: **90W**
Protection: **IP X6**

JAN	FEB	MAR	APR	MAY	JUN	JUL	AUG	SEP	OCT	NOV	DEC
2022	2023	2024	2025	2026	2027						

- Do not throw away. Send to a recycling centre.
- The packaging should be disposed of as reusable waste.
- Metals should be understood as old material.
- Plastic products should be treated as recycling material.
- Waste management should be as established in each country.
- Please ask the local administration about recycling and rubbish.

Read the instruction manual before use.

Maximum load 250 Kg.

Safety warning and instructions.

1.2. Intended product use

The SILVER + (TG-1200) stretcher has a motorised design for the transport and support of people in hospital and pre-hospital environments. The battery-powered electro-hydraulic lifting system has been developed to reduce the effort required by the caregiver when lifting and lowering the stretcher. This stretcher has a maximum capacity of 250 kg and requires trained professional operators to use it, including emergency medical services, care home medical staff and first aid medical teams.

The SILVER+ stretcher is designed for patient transport and is not intended for prolonged stays or for use as a hospital bed, nor for use in devices that modify atmospheric pressure, such as hyperbaric chambers.

1.3. Liability and warranty

The stretcher should be checked when handed over to the relief organisation. All of the functions should be explained in detail. The relief organisation is responsible for training all employees in the correct use of the equipment.

The product warranty is 24 months from the date indicated on the purchase invoice. (See Warranty Conditions and Warranty Certificate enclosed with the product).

The Warranty will not cover breakdown due to incorrect installation, mishandling or misuse of the product. Repairs must be carried out by a service technician authorised by Kartsana or its authorised representative.

The manufacturer shall not be liable for any damage to the stretcher caused by the use of non-original Kartsana products.

Changes or modifications made to the unit that have not been expressly approved by Kartsana may result in errors in the equipment's operation.

1.4 Specifications

Maximum workload (unassisted)	250 kg
Main regulations	UNE-EN-1865-2 EN-1789+A1 UNE-EN 60601-1 / UNE-EN 61000-6
Articulation/tilting of the backrest	From 0° to 80°
Articulation/tilting of the leg	Simple articulation
Diameter of rotating wheels	160 mm
Total length ¹	1965 mm
Total width	568 mm

Minimum and maximum height ²	350 mm - 1141 mm
Loading height ³	Max. 791 mm
Recommended operators for loading/unloading an occupied stretcher	1 operator or 2 operators in the event that the stretcher is used by a very large person
Compatible restraining systems	Rail Silver+ (R-1200)

¹ The overall length can be shortened to 1607 mm (with the headrest extension retracted to a minimum and the backrest folded to 80°C - position used in tight manoeuvres, lifts, etc.).

² Height from the floor surface to the structural profiles of the patient area of the stretcher at its minimum and maximum height position.

³ For patients weighing more than 200Kg, it is recommended to perform the movements with the stretcher at a low working height.

Kartsana reserves the right to change the specifications without prior notice.

1.5 Attention

Throughout this manual, there are 3 types of symbols. As a **warning symbol**, we have ISO 7010-W001. The **prohibition** is indicated by ISO 7010-P001 and the **obligation** by ISO 7010-M001.

Using these symbols, we provide information on important safety measures for the correct use of the stretcher and to avoid possible accidents.

Warnings alert the reader to situations that, if not avoided, may lead to death or serious injury. Furthermore, prohibitions indicate actions that must not be carried out under any circumstance, and obligations indicate measures and/or actions that must be taken to ensure the proper use of the product and to maintain safe functionality.

1.6 Summary of safety precautions

Read and follow the warnings and cautions on these pages carefully. Maintenance and repairs may only be carried out by qualified personnel.

NOTICES

- Improper use of the stretcher may cause injury to the patient or the handler. Use the stretcher only as described in this manual.

- Do not modify the stretcher or any of its components. Modification of the device may result in an unpredictable operation that could cause injury to the patient or operator. Modification of the device voids the warranty.

- The emergency vehicle where this stretcher is used, must have a compatible restraining system installed.

- The emergency vehicle fastening system compatible with our stretcher must be correctly installed for correct attachment of the stretcher. Otherwise injury to patient, operator or vehicle could occur.

- Do not attempt to operate the stretcher when it is secured to its fastening system in the vehicle.

- Request the correct installation of the fixing rail to the vehicle floor. The position of this rail will determine the correct functioning of the stretcher. After installation, check that the wheels and forks of the stretcher are raised and do not come into contact with the bumper of the vehicle.

- To avoid risk of electric shock, do not attempt to open the battery for any reason. If the battery case is cracked or damaged in any way, do not insert the battery into the

charger. Please return the batteries to be replaced with the correct ones.

- Do not remove the battery when the stretcher is charging inside the ambulance.
- Avoid direct contact with the battery or battery casing when wet, as this may cause injury to the patient or operator.
- Check batteries and chargers for damage before each use.
- Practice changing the height and load positions of the stretcher until you have fully mastered the operation of the product. Improper use can lead to injury.
- Do not allow untrained assistants to help you use the stretcher. Inadequately trained technicians and assistants can cause injury to the patient and to themselves.
- Do not climb on the base of the stretcher as this may damage the product, unbalance the patient and cause injury.
- Transporting the stretcher on its side could cause it to tip over and result in damage to the product and injury to the patient. Transporting the stretcher in a lowered position reduces the chances of the stretcher tipping over.
- Incorrect fastening of the stretcher may result in injury. Avoid putting hands, fingers or feet on any moving parts of the stretcher. To avoid injury, when raising and lowering the stretcher, be very careful when placing hands and feet close to the moving structure (the structural scissors).
- Always use all the restraining belts to secure the patient on the stretcher. If not properly restrained, the patient could fall off the stretcher and injure him or herself.
- Never leave a patient unattended on the stretcher, as they may become injured. When transporting a patient, hold the stretcher firmly in place.
- Do not lock the lower wheels if you want to move the patient. The stretcher could tip over if an attempt is made to move the stretcher with the brakes locking the wheels engaged, with the risk of injury to the patient or operator, and damage to the stretcher.
- Side rails are not designed to be used as a means of restraint.
- Due to the electro-hydraulic raising and lowering of the stretcher, the equipment used for patient monitoring may be temporarily affected. For best results, patient monitoring should be done when the stretcher is not moving.
- Architectural barriers such as curbs, steps or uneven terrain can cause the stretcher to tip over, potentially injuring the patient or operator.
- The SILVER+ (TG-1200) stretcher is designed to be compatible with the R-1200 rail, and it is the responsibility of the operator that these products work together.
- If the stretcher is occupied by a person of large size, we recommend two operators to handle the stretcher and secure the patient.

- The higher the operator has to lift the stretcher, the more difficult it will be to support the weight. If the operator is very short or if the patient is too heavy to be lifted safely, assistance may be required to load the stretcher. A short operator will have to raise his or her arms higher in order to unfold the lower part of the stretcher.

- Do not install or use the stretcher brakes with excessively worn wheels. Installation or use of the brakes on a wheel with a diameter of less than 160 mm could compromise the holding capacity of the brakes, with the risk of injury to the patient or operator and damage to the stretcher or other equipment.

- When cleaning, use appropriate personal protective equipment (protective goggles or respirators) to avoid the risk of inhaling infectious organisms.

- Some cleaning products are corrosive and can damage the product if used incorrectly. If the products described above are used to clean Kartsana equipment, steps should be taken to ensure that the stretchers are cleaned with clean water and thoroughly dried after cleaning. Failure to properly rinse and dry the units may leave a corrosive residue on the surface of the units, which may cause premature corrosion of critical components.

- Failure to properly clean or improperly dispose of a contaminated mattress or other components increases the risk of exposure to blood-borne pathogens and could result in injury to the patient or operator.

- Leaking pressurised fluids can penetrate the skin and cause serious injury. To avoid this danger, release the pressure before disconnecting hydraulic or other lines. Before applying pressure, tighten all connections. If an accident occurs, consult a doctor immediately. Any fluid that has penetrated the skin must be surgically removed within a few hours, otherwise it may become gangrenous. Doctors unfamiliar with this type of injury should refer the patient to a specialised medical centre.

- To avoid the risk of injury, wear gloves when checking hydraulic connections for leaks.

- Take special precautions for electromagnetic compatibility (EMC) when using medical electrical equipment such as the Silver+ TG-1200 stretcher. The use of portable and mobile radio frequency (RF) communications equipment may affect the operation of the Silver+ TG-1200.

- The use of accessories, transducers or cables other than those specified, with the exception of transducers and cables marketed by Kartsana as replacement parts for internal components, may result in increased emissions or decreased immunity of the Silver+ TG-1200 stretcher.

PRECAUTIONS

- Changes or modifications to the unit not expressly approved by Kartsana may void the user's authority to operate the system.
 - The stretcher can be adjusted to any loading height position. Before the stretcher is put into service, establish the required loading height of the stretcher.
 - Adjust the loading height of the stretcher to the appropriate standing height before use.
 - Installation of the stretcher-compatible fixing rail must be carried out by a certified mechanic familiar with the ambulance structure. Consult the vehicle manufacturer before installing the rail and ensure that the installation does not damage or interfere with the vehicle's brake, oxygen or fuel hoses, fuel tank or electrical wiring.
 - Use the battery and charger only as described in the user/maintenance manual.
 - The stretcher cannot be used with an AC adapter.
 - If charging the battery in the ambulance, place the charger in a locked cabinet out of reach of the patient during transport.
 - Check that the battery is fully charged before putting it into service. A discharged or flat battery may result in poor stretcher performance.
 - Before putting the stretcher into operation, remove all obstacles that may interfere and cause injury to the operator or patient.
 - When unloading the stretcher from the patient compartment, make sure that the rotating wheels are placed securely on the floor, otherwise the product may be damaged.
 - Remove the battery if the stretcher is not to be used for an extended period of time (more than 24 hours).
 - The wheel brakes are uniquely designed to help prevent the stretcher from rolling away when unattended and to make it easier to transport patients. The wheel brakes may not offer sufficient resistance on all surfaces or with weights of the wheels..
 - Ensure that the safety belts do not become entangled in the base frame when raising and lowering the stretcher.
 - Do not store any objects under the stretcher mattress. Storing objects under the mattress may interfere with the functioning of the stretcher.
 - Do not use steam or ultrasonic cleaners to clean the unit.
 - The maximum temperature of the water must not exceed 80°C.
 - Allow the stretcher to air dry.
 - Dry all rotating wheels and interface points with a towel.
- Failure to comply with these instructions may result in some or all of the warranties becoming void.
 - Before cleaning the stretcher, remove the battery.
 - A preventive maintenance programme should be established for all Kartsana equipment. Depending on the product's frequency of use, preventive maintenance may need to be carried out more frequently. Close attention must be paid to the safety features, which include, among others: the hydraulic mechanism and electrical controls.
 - For further information on maintenance, please refer to the information on preventive maintenance.
 - Improper maintenance may cause injury or damage to the product. The stretcher should be maintained as described in this manual. Follow the maintenance procedures and use only Kartsana approved parts. The use of unapproved parts or procedures may result in an unpredictable operation and may cause injury and void the warranty.
 - Failure to use authorised spare parts, lubricants, etc. may damage the stretcher and void the product warranty.
 - Hydraulic lines, hoses and fittings can fail or become detached due to physical damage, kinks, wear and environmental exposure. Regularly check hoses and lines to prevent damage to the stretcher. Tighten any loose connections.

NOTES

- Loose elements or dirt accumulation on the floor of the patient compartment can interfere with the operation of the rail and the stretcher restraining system. Keep the floor of the patient compartment clean.
- Kartsana works constantly to improve the design and quality of its products. Therefore, although this manual contains the most current information available, there may be minor discrepancies between the stretcher and this manual. If you have any questions, please contact Kartsana's customer service or technical service area.
- Kartsana recommends that, prior to installation, a certified mechanic should plan the placement of the mounting rail on the inside of the emergency vehicle.
- Automatic charging is only possible with batteries supplied by Kartsana.
- With the Silver+ (TG-1200), use only Kartsana approved batteries.
- If the batteries are not on the charger, they will slowly lose power.

2. INSTRUCTION MANUAL

2.1. Technical characteristics of the SILVER+ (TG-1200) stretcher



Do not damage the areas where the stretcher mechanisms are located in order to avoid malfunctioning of the stretcher mechanisms.


Stretcher weight: 73 Kg

Dimensions and measurements in mm


Stretcher with folded legs:

2.2. Vehicle configuration

These installation instructions are indicated for the attachment of the SILVER+ (TG-1200) set with the Kartsana SILVER+ (R-1200) rail.

The correct use and operation of the SILVER+ (TG-1200) stretcher depends on the correct installation of the fastening system (rail R-1200).

The fixed part of the rail shall be fixed at a **maximum distance of 50 mm (A)** from the rear door of the ambulance.

Careful consideration should be given to the distance/type of **bumper** that the ambulance model is fitted with, so that if the bumper protrudes more than 180 mm it will have to be adapted to allow the stretcher legs to be raised correctly.

Depending on the ambulance model and the type of rear bumper, the fixing dimensions may vary and/or be adapted.



The height of the ambulance floor in relation to the ground **must not exceed 70 cm.**

2.3 Operation and handling of the stretcher

2.3.1 Fitting the batteries

The first step in working with the Silver+ is to fit the 28 V battery that powers the system. It is positioned at the rear (at the patient's feet) in the operator's monitoring area.

To remove the battery, press the red buttons simultaneously, and pull to unclip the battery from its connector.

The battery's autonomy in no-load operation (no load) should be about 190 cycles.

Information note: The battery is delivered with a minimum charge. Before operating the stretcher, the batteries must be charged with the supplied charger. After that, we will be able to work normally.


If the battery becomes excessively discharged, the battery must be charged with the external battery charger. Do not attempt to charge the excessively discharged battery with its Rail Silver+ charging point (R-1200).

To charge the batteries outside the ambulance, the specific charger is included with the stretcher.

2.3.2 Setting up and control of the stretcher

To avoid possible injury to the patient and attendants in the event of an accident, do not place objects in the path of the moving parts of the stretcher. It is also recommended to avoid protruding elements inside the ambulance.

The system starts up once the battery is connected to the stretcher and the power switch is turned on. The number of LEDs lit on the control unit will indicate the charging status of the control unit.

If an error occurs during the self-control process, this will be indicated by the error LEDs (see error table in chapter 7).

Functions of the buttons on the control panel:

The buttons on the keypad function as follows:


The stretcher goes into low-power Stand-By mode if no key is pressed for 30 minutes. Simply press any button to turn the stretcher back on.

2.3.3. Adjustable head

Operate the red lever under the head and move the head to the desired position.



Then release the lever so that the head is locked in the chosen position.

2.3.4. Adjustable leg

Button to adjust the inclination of the legs



By pressing this button and moving the leg rest manually, it is possible to adjust the positioning of the legs to the desired inclination. Once the inclination is correct and the button is released, the clamp is automatically locked in the required position.


2.3.5. Safety railing

The safety rail has a red handle that locks the rail. To fold down the railing in the longitudinal direction, just press this handle as shown in the figure below. To re-lock it in the safety position, it must be moved manually to this position and it locks automatically.

This would allow the stretcher to be left with the stretcher rails folded lengthwise:

The railing can also be folded in the outward direction. To do this, pull both handles at the same time and follow the movement of the handrail outwards.



2.3.6. Drip holder

The stretcher is equipped with a drip holder on the right side profile. This is screwed in place with 2 screws to the stretcher frame profile.

The length of the drip tube can be lengthened or shortened by pressing the red lever at the top of the black cylinder.

It can be positioned on both sides of the stretcher (right and left).


Folding

Dimensions:
467 mm x 60 mm
(extended is **390 mm longer**)

Weight:
0.68 kg


Maximum permissible load: **6 Kg**

Extended


2.3.7. Safety belts

Head belts:

Place the fasteners (C) in the position shown in the figure above. Pass the male buckle (D) between the latches (C) until it locks into buckle B. Once the assembly is secured, centre and tighten the assembly by adjusting the belt tension at the ends A, E, F and/or G. (The latches C do not appear on the foot belt).

Lap belts:

Place female buckle A and male buckle B in the position shown in the side figure and insert the latter into the groove of the former. Once the assembly is secured, centre and tension the assembly by adjusting the tension of the belt using end C.


It is recommended that the patient is always secured with safety belts during the time he or she is on the stretcher.

Belt knot

2.3.8. Brake lever for rear wheels


By **pressing** the lever indicated in the drawing with your foot, the rear wheels are locked by the brake. These assist in keeping the stretcher static, they should never be used as dynamic brakes.

2.3.9. Lever for locking the front wheels

Pressing the **red** central lever allows the wheels to turn freely.

Pressing the **green** rear lever blocks the free rotation of the wheels.

The front wheels must always be locked when loading the stretcher onto the rail. Failure to do so could damage important parts of the stretcher.

Do not operate the brake lever if there is a patient on the stretcher. The stretcher may tip over if moved with the brakes applied, with the risk of injury to the patient or operator, and damage to the stretcher.

Never leave a patient unattended on the stretcher, as this could result in injury; hold the stretcher securely when a patient is on it.

Do not install or use the stretcher brakes with excessively worn wheels. Installation or use of brakes on wheels with a diameter of less than 150 mm could result in malfunction and reduced braking power, with the risk of injury to the patient or operator and damage to the stretcher.


The wheel brakes are designed only to help limit the movement of the stretcher when it is not to be moved and to facilitate the transport of patients. The wheel brakes may not offer sufficient resistance on all surfaces or with weights of the wheels..

For a patient weighing more than 200 kg, it is recommended that a low position be used when transferring the patient, and the stretcher can be raised to the required position when the patient is transferred into the ambulance.

2.3.10. Mounting of the SILVER+ stretcher with the R-1200 rail

The maximum range of the stretcher should only be adjusted **the first time** the stretcher is loaded onto the rail. If moving to another ambulance of a different height, **it will have to be readjusted**.

To set the maximum height of the stretcher's upward trajectory, the stretcher should be positioned in front of the deployed rail on the outside of the ambulance. Position the front wheels half a centimetre above the profile where the small load wheels slide. (*see side picture*).


- 1** The stretcher must be kept in this position when removed from the rail.
- 2** Then look for the adjustment cam on the top of the stretcher near the shock absorber attachment on the headrest.
- 3** Rotate the cam until the end-of-range light comes on.
- 4** Once in position, tighten the screw 8 Nm to fix the cam in its final position.

Thereafter, the stretcher shall be raised to the maximum level. Check the loading manoeuvre and readjust if necessary.


On uneven terrain or changes in gradient, the stretcher can be raised to its maximum range by pressing the right red and yellow button at the same time.

Loading the stretcher on the rail



1 Face the stretcher, centring the front wheels with the rail castings.

2 Press the left button to raise the legs. Keep pressing it down until it no longer rises.

3 Press the blue lever to release the pin lock.

4 Slide the stretcher along the rail guides until the front wheels stop and the third point engages.

To ensure that the stretcher-rail assembly is securely attached, make sure that the release lever is in the horizontal position. The 3rd point being inside its housing.

If the stretcher receives power from the ambulance and the stretcher is in the correct position, battery charging starts. During the charging process, the battery LEDs will flash in a sweeping pattern, indicating that charging is in progress.


Once charging is complete, LED A turns on and the four red battery LEDs turn off.


Unloading of the stretcher when on the rail

To unload the stretcher, the procedure is as follows:


- 1 Pull the lever upwards to extend the rail forwards.

- 2 Push the red lever on the rail downwards.


- 3 Slide the stretcher until it stops.

- 4** Unfold the legs by pressing the right button on the control unit.

- 5** Pull the blue lever on the right to unlock the mechanism and pull out.

2.3.11. Transferring the patient to the stretcher

Transferring the patient to the stretcher

1. Move the stretcher towards the patient.
2. Place the stretcher next to the patient and raise or lower it to the patient's level.
3. Lower the side rails and open the restraining belts.
4. Transfer the patient to the stretcher using standard emergency medical services procedures.
5. Use all the restraining belts to secure the patient on the stretcher.
6. Adjust the backrest and footrest as required.

Moving the stretcher:

1. Ensure that all the restraining belts are secured above the patient.
2. There should be one operator at the lower end and one at the upper end of the stretcher at all times when moving the stretcher with a patient.
3. Approach door thresholds or other low barriers directly and lift the wheels over the obstacle separately, never all four at once.

Appropriate lifting techniques:

When lifting the stretcher and the patient, follow these five basic guidelines to help avoid injuries:

- ◇ Keep hands close to the body.
- ◇ Keep the back straight.
- ◇ Coordinate movements with your partner and lift with your legs.
- ◇ Avoid turning around.
- ◇ Always use the stretcher as described in this manual.

3. FAULTS

In the event of a hydraulic system failure, pull the red handle located next to the batteries at the rear of the vehicle. The stretcher will lower to the lower limit and remain in this position. Please note that if the stretcher is fully extended to its mechanical limit, manual lowering using the trigger will **NOT** work.

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USE ONLY IN CASE OF HYDRAULIC OR ELECTRIC FAILURE

Please use cleaning agents according to the manufacturer's instructions. As shown in the diagram, the chassis consists of 2 folding handles.

It is recommended that the emergency lift is carried out by 2 people, one on each side.

After performing the normal loading manoeuvre as described above, the stretcher should be pushed into the rail until it reaches the stop.

4. CARE AND MAINTENANCE OF THE PRODUCT

4.1. Contact information

If you have any questions regarding the use, maintenance or installation of a KARTSANA product, please contact our Technical Service Area (SAT). Any serious incident involving the device must be reported to the manufacturer and to the competent authority of the Member State where the user and/or the patient are established.

Telephone: +34 93 715 86 72
Email: sat@kartsana.com
Address: Narcís Monturiol, 34
08192 Sant Quirze del Vallès,
Barcelona (Spain)

4.2. Environmental conditions

Both for use and storage:

Recommended temperature: from -10°C to +60°C

Relative humidity: 5 – 95%

4.3. Cleaning

To clean the stretcher, use a cloth, rag or similar to wipe the surface to be cleaned with products that do not damage the surface of the materials. If the surface to be cleaned is already lubricated, re-grease it after cleaning has been completed.

Please use cleaning agents according to the manufacturer's instructions.

Do not wash the stretcher with pressurised water, take special care with electrical elements, housing, electronic board, keyboard, LEDs and motor.

Remove the batteries before washing and dry the terminal contacts afterwards.

The manufacturer shall not be liable for any damage or anomaly caused by a cleaning product that may damage the surface of the stretcher materials.

For reasons of hygiene and preservation of the product components, the manufacturer recommends cleaning after each use.

4.4. Maintenance

4.4.1 Precautionary maintenance

Correct and regular maintenance guarantees the durability of the product.

Creating a maintenance plan, including regular reviews and establishing a responsible person to carry them out, is recommended.

The person in charge must adhere to the following **requirements**:

- Possess technical knowledge related to the product and the maintenance procedures described in these instructions.
- To have the support of qualified, trained and knowledgeable staff who are familiar with the operation of the product.
- Use original or manufacturer-approved components, spare parts and accessories.
- In accordance with the instructions of the European Council Directive 93/42/EEC which establishes the obligation of the purchaser to provide on request the above mentioned after-sales care record for the purpose of product traceability, records of all product maintenance operations shall be kept.

Check for proper operation **before use**. Prior to service, the following should be checked:

- Functionality.
- Fastenings and hardware.
- Parts that move, such as wheels, belts and mattresses.
- Sensitivity of springs.
- Electrical functions: Full range of motion, lighting, etc.

If the device does not appear to be suitable for correct and safe use, **it should be taken out of service** until the device is repaired or restored to full working order.

Do not modify the structure of the device for any reason, as this may cause serious harm to patients and/or operators.

Lubrication intervals of moving parts:

- 1 – 30 services per month: every 3 months.
- 30 – 50 services per month: every 2 months.
- More than 150 services per month: every month.

All maintenance other than lubrication, tightening of nuts and bolts and ordinary cleaning must be carried out by KARTSANA or an authorised service centre.

4.4.2. Service maintenance

The person entrusted with the product review shall ensure the following basic requirements:

- ◇ Adequate knowledge of the product, its technical/constructional characteristics, control points and final tests, packaging, preservation and handling.
- ◇ Knowledge of all product functions and any possible risks, malfunctions or failures of the product.
- ◇ Possession of all necessary tools to carry out any technical service or minor repairs.
- ◇ Possession (or ability to procure) original manufacturer's or manufacturer's authorised spare parts.
- ◇ Use or support of specialised technical personnel trained by the manufacturer to service the product in question.
- ◇ Keeping a record of all maintenance operations carried out on the product, in accordance with the instructions of the European Council Directive 93/42/EEC which establishes the obligation of the purchaser to provide on request the above mentioned after-sales care record for the purpose of product traceability.

4.5. Transport and storage

Before transporting the device, ensure that it is properly packaged and secured to prevent damage during transportation.

Keep the original packaging for possible further transport. Any damage caused to the device during transport is not covered by the warranty.

Repairs and replacements of damaged parts are the responsibility of the customer. The device should be stored in a dry place.

When storing, do not place any heavy objects on top of the product. It must not be used as a support for any kind of object.

4.6. Disposal

When the product is no longer in a usable condition, do not throw it away but develop a take-back strategy for recycling by an authorised dealer. The packaging shall be managed as reusable waste. Metals are fully recyclable. Plastic products should be treated as recycling material. Waste management should be as established in each country. If you have any questions, please ask your municipal administration.

5. SPARE PARTS

To consult the spare parts for this product, you can visit the **following link:**

https://www.kartsana.com/global/Spare_parts_TG1200.pdf

If you have any questions, please contact our technical service.

6. STRETCHER ACCESSORIES

C-090 Drip tube

C-006 Fixed head extension

C-014 Folding head extension

C-040 Corpuls support


C-101 Instrument table

07-0112 Battery charger

C-095 Oxygen bottle holder

P-370 Mattress

07-0110 Battery

7. SYSTEM ERRORS

When an error occurs, the **stretcher becomes inactive** showing:

- Folded stretcher indicator flashing.
- The marking of the error code by a combination of the four yellow error LEDs.

To exit the Error status:
Press the red UP button + the red DOWN button + the yellow LIGHT button
at the same time for 2 sec.

7.1. Error codes

ERROR NAME	CODE	LEDs			
		A	B	C	D
Battery Charging Error due to NOT exceeding minimum voltage	1			●	
Charge Intensity Error VERY HIGH Batteries	2			●	●
Battery Charging Error due to lack of Ambulance Voltage	3		●		
DISCHARGED or VERY LOW Batteries	4		●		●
Ambulance Switch Error	5		●	●	
Charger Error due to lack of battery, or stretcher disconnected	6		●	●	●
Load Error due to exceeding Duty max in PWM	7	●			
Error in Upward Movement due to excessive time	8	●			●
Error in Downward Movement due to excessive time	9	●		●	
Error in Motor due to exceeding maximum current	10	●		●	●
Error in Motor due to low voltage	11	●	●		

7.2. Troubleshooting

The following table summarises the possible errors and their respective solutions.

In addition, it is advisable **to remove the battery if the stretcher is not going to be used** for long periods of time.

TROUBLESHOOTING			
CODE	DESCRIPTION OF THE ERROR	POSSIBLE CAUSE OF THE ERROR	ERROR SOLUTION
1	The battery voltage is less than 3V. The battery is disconnected or damaged.	Battery disconnected	Connect the batteries to the battery connector correctly.
			Check that the battery cables are properly connected to the electronics board.
		Damaged battery	Replace with a new battery.
2	The battery voltage is higher than 34V.	An overload has occurred in the battery recharging process. (Use only the charger supplied by KARTSANA or the stretcher for battery recharging).	Leave the battery idle, preferably in a cool place, for 3 to 4 hours (or until it does not show error 2 when connected to the stretcher) to allow it to self-discharge. Check that the batteries have not been damaged internally by discharging them by moving the stretcher up and down. It should provide at least 15 cycles before discharging.
		Damaged battery	Replace with a new battery.
3	The stretcher does not receive voltage from the ambulance.	The stretcher-rail plug is not connecting correctly	Check that the two contacts of the stretcher plug and the rail are touching each other and that there is no foreign body between them.
		Ambulance battery is discharged or damaged.	Charge the battery or replace it with a new one, respectively.
		Power cables incorrectly connected either on the ambulance, rail or stretcher.	Check connection of the cables. Recommended order of testing: ambulance, rail and stretcher.
		One of the copper contacts of the stretcher plug has not come out of its housing completely in the plastic part.	Check that the plug contacts move smoothly in and out of the socket.
4	The battery voltage is less than 18V.	Discharged battery	Charge the battery in the charger and not in the vehicle.
5	The stretcher is receiving power from the socket, but the position sensor has not been activated.	The stretcher has not fully entered the rail.	Position it correctly inside the rail.
			Adjust the position sensor.
6	Battery is not detected.	Disconnected batteries	Push them fully into the clamping plate until you hear a clicking sound.
		Damaged battery	Replace with a new battery.
		Very low battery	Charge with external charger.

7	Failure in charging process.	Contact technical support	Contact technical support.
8	More than 9 seconds have passed in an upward motion.	Discharged battery	Charge the battery
		Hydraulic system failure.	Contact technical support.
		The mechanism does not move freely.	Check for excessive friction or interference between the mechanical parts.
9	More than 9 seconds have elapsed in downward motion.	Discharged battery	Charge the battery
		Hydraulic system failure	Check for possible oil leaks. Contact technical support.
		The mechanism does not move freely	Check for excessive friction or interference between the mechanical parts.
10	Current greater than 48A.	Overloaded stretcher	Do not load the stretcher with loads of more than 250 Kg.
		Discharged battery	Charge the battery
		The mechanism does not move freely.	Check for excessive friction, blockage or interference between the mechanical parts.
11	Lower voltage of 15V.	Overloaded stretcher	Do not load the stretcher with loads of more than 250 Kg.
		Discharged battery	Charge the battery
		The mechanism does not move freely.	Check for excessive friction, blockage or interference between the mechanical parts.

7.3. Motor protection

The system switches off in the following situations:

- By Time-out of 12 sec (*error 8 and 9*)
- For Intensity greater than 60A for more than 0.5 sec (*error 10*)
- For Motor Battery Voltage less than 18 Volt (*error 4*)

7.4. Ambulance entry

If there is ambulance entry, 4 different cases can occur:

- 1) Stretcher on without ambulance voltage:** the LEDs will be off waiting for vehicle voltage. After 10 seconds the LEDs on the stretcher flash (*error 3*). They will stop flashing once the stretcher receives ambulance voltage.
- 2) Stretcher off without ambulance voltage:** No error appears, but the stretcher is not charging. Ideal scenario for prolonged breaks without using the stretcher.
- 3) Stretcher off in ambulance:** LEDs off waiting for battery detection. After 3 seconds the LEDs on the stretcher flash (*error 6*). They will stop flashing when the stretcher is turned on.

4) Stretcher switched on in ambulance: If the stretcher receives power from the ambulance and the sensor position is correct, the battery charging process starts. During charging, the battery LEDs will flash in a sweeping pattern, indicating the operation in progress. Once charging is complete, LED A lights up and the four red battery LEDs turn off.

7.5. Battery charging

When entering the ambulance, the battery is tested for 10 seconds to check that its voltage is between 24 and 34 Volts. Furthermore, 100 mA current is applied to check that the battery is connected. The battery charging process consists of 2 phases.

For the first phase, one of the following conditions must be met:

- The previous upload was incomplete
- The battery is discharged
- The battery has a voltage of less than 25.6 Volts.

PHASE 1

It consists of 3 consecutive steps (must be completed to determine that the battery is charged):

1) Start of charging:

Charge with a steady current of 750mA until one of the following conditions is met:

Voltage greater than 24 Volts Maximum time 20 minutes Ambulance departure Error by:	Duty Max Voltage greater than 30.4 Volts Current greater than 6.0A Does not exceed the minimum voltage of 17.6 V within 30 minutes
--	---

2) Fast charging:

Charging at a constant current of 1.8A until one of the following conditions is met:

The voltage of 28.8 Volts is reached Maximum time 120 min Ambulance departure Error by:	Duty Max Voltage greater than 30.4 Volts
--	---

3) Full charge:

Charge at a constant voltage of 28.8 volts until one of the following conditions is met:

Maximum time 60 minutes
Ambulance departure
Error by:

Duty Max
Voltage greater than 30.4 Volts
Current greater than 6.0A

For the second phase, it is only required **that no errors have occurred.**

PHASE 2.

Consists of 1 step, in this case, the maintenance charge does not need to be completed:

1) Maintenance charge:

Charge with a constant current of 80mA until one of the following conditions is met:

Maximum time 10 minutes
The saturation voltage of 29.04 Volts is reached.
Ambulance departure
Error by:

Duty Max
Voltage greater than 30.4 Volts

If the battery falls below 20 Volts, the battery must be charged by the battery charger and not by ambulance (*error 4*)

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Model: **TG-1200**
 Commercial name: **SILVER+**
 Voltage: **24V DC**
 Max. Power Input: **90W**
 Protection: **IP X6**

	JAN	FEB	MAR	APR	MAY	JUN	JUL	AUG	SEP	OCT	NOV	DEC
	2022	2023	2024	2025	2026	2027						

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