Serie

JÚPITER

CARRO-CAMILLA

TG-880 2 / 4
MTR-120 / R2

KARTSANA

Life-moving innovation

C/ Narcís Monturiol, 34
08192 Sant Quirze del Vallès
Barcelona (España)
+34 93 715 86 72
info@kartsana.com
www.kartsana.com

Manual de instrucciones

**Camilla mecánica TG-880 2 / 4
MTR-120 / R2 / MT-118IN**

Instruction manual

**Mechanical stretcher TG-880 2 / 4
MTR-120 / R2 / MT-118IN**

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**Brancard mécanique TG-880 2 / 4
MTR-120 / R2 / MT-118IN**

ESPAÑOL

ENGLISH

FRANÇAIS

JÚPITER CARRO-CAMILLA

Camilla mecánica

Modelo TG-8802 / 4 | MTR-120 / R2

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
Mechanical stretcher

Model TG-8802 / 4 MTR-120 / R2

Please read these maintenance instructions carefully before using the product and keep them for future reference.

CONTENTS

1. INTRODUCTION	34	3. CARE AND MAINTENANCE OF THE PRODUCT	57
1.1 Preface	34	3.1 Contact details	57
1.2. Intended product use	34	3.2 Environmental conditions	57
1.3. Liability and warranty	34	3.3 Cleaning	57
1.4 Specifications	35	3.4. Maintenance	58
1.5 Symbols	36	3.4.1 Precautionary maintenance	58
1.6 Summary of safety precautions	36	3.4.2. Service maintenance	59
2. INSTRUCTION MANUAL	38	3.5 Transport and storage	59
2.1. Technical specifications and dimensions of the TG-880 + MTR stretcher-trolley	38	3.6 Disposal	59
2.1.1. Stretcher-to-trolley locking mechanism	40	4. SPARE PARTS	60
2.2. Technical specifications and dimensions of the MTR-120 series stretcher	41	5. STRETCHER ACCESSORIES	61
2.2.1. Operating the stretcher	43		
2.2.2. Removable handles	43		
2.2.3. Adjustable backrest section	44		
2.2.4. Adjustable leg section	44		
2.2.5. Safety railing	45		
2.2.6. IV pole	46		
2.2.7. Safety belts	46		
2.2.8. Removing and fitting the stretcher belts	47		
2.3. Incubator version	48		
2.3.1. Product description and dimensions	48		
2.4. Technical specifications and dimensions of the TG-880 series trolley	50		
2.4.1. Rear controls for locking/unlocking the rear wheels	52		
2.4.2. Controls for operating the stretcher legs	52		
2.4.3. Intermediate positions of the stretcher trolley	53		
2.4.4. Brake lever	55		
2.4.5. Mounting the trolley on the ambulance stretcher holder	55		

1. INTRODUCTION

1.1 Preface

The **Jupiter series of stretcher trolleys (TG-880 + MTR-120)** have been specifically designed for patient rescue and transport purposes. The product complies with the MDR (EU) Regulation 2017/745. It has also been tested in accordance with the UNE-EN ISO 1865-3 and UNE-EN 1789 standards.

This series comprises two different trolley models and two different stretcher models. When referring to the trolleys, we will use the designation **TG-880**. When referring to the stretchers, we will use the designation **MTR**. This series of trolley-stretchers must be secured within the ambulance using the range of rails specified in this manual.

All information on treatment, disinfection and maintenance is based on our current experience and understanding.

We reserve the right to make variations to the technical specifications of the stretcher in order to improve it.

1.2. Intended product use

This **TG-880 + MTR** set functions as a motorised stretcher designed for the transport and support of people in hospital and pre-hospital settings. It comprises two separate products which, although they function together as a single unit, still have independent functions. Its reversible function allows for greater versatility in patient care, both inside and outside the ambulance. The connection between the two products is convenient and secure.

The mechanical lifting system, with 8 intermediate positions, is designed to adapt to the required working conditions. This stretcher has a maximum capacity of 280 kg and requires trained professional operators to use it, including emergency medical services, care home medical staff and first aid medical teams.

The **TG-880 + MTR** stretcher trolley is designed for patient transport and is not intended for prolonged stays or for use as a hospital bed, nor for use in devices that modify atmospheric pressure, such as hyperbaric chambers.

1.3. Liability and warranty

The product must be checked at the time of delivery to the emergency response organisation. All of the features must be explained in detail. The emergency response is responsible for training all employees in the correct use of the equipment.

The product has a warranty period of 24 months from the date indicated on the purchase invoice. (see *Warranty Conditions* and *Warranty Certificate* enclosed with the product).

The warranty will not cover faults resulting from improper installation, misuse or mishandling of the product. Repairs must be carried out by a service technician authorised by Kartsana or its authorised representative.

The manufacturer will not be liable for any damage to the stretcher caused by the use of non-original Kartsana products.

Changes or modifications made to the unit that have not been expressly approved by Kartsana may result in errors in the equipment's operation.

1.4 Specifications

SPECIFICATIONS	TG-880 2	TG-880 4	MTR-120	MTR-120 R2
Maximum workload	280 kg			
Main regulations	UNE-EN-1865-3 EN-1789 EN-1789+A1			
H1 Height (mm) ¹				
H2 Drive wheel height (mm)				
H3 Folded height				
W Width (mm)				
L Length (mm)				
Weight (kg)				
Recommended operators for loading/unloading an occupied stretcher	1 operator or 2 operators if the stretcher is being used by a patient requiring assisted handling due to body weight.			
Compatible securing systems	R-419L R-450N R-450N INOX R-800L			
Front wheels	Fixed	Lockable swivel castors	-	
Leg section	-		Double articulated	Single

¹For patients weighing more than 200 kg, movements should be carried out with the stretcher set to a low working height. Kartsana reserves the right to change the specifications without prior notice.

1.5 Symbols

Throughout this manual, there are 3 types of symbols. As a **warning symbol**, we refer to ISO 7010-W001. The **prohibition** symbol is indicated by ISO 7010-P001 and the **obligation** by ISO 7010-M001.

We use these symbols to provide information on important safety measures for the correct use of the stretcher and to prevent accidents.

Warnings alert the reader to situations which, if not avoided, may lead to death or serious injury. Furthermore, prohibitions indicate actions that must not be carried out under any circumstance, and obligations indicate measures and/or actions that must be taken to ensure the proper use of the product and to maintain safe functionality.

1.6 Summary of safety precautions

Carefully read and follow the warnings and precautions set out in these pages. Maintenance and repairs may only be carried out by qualified personnel.

NOTICES

- Improper use of the stretcher may cause injury to the patient or the handler. Use the stretcher only as described in this manual.
- Do not modify the stretcher or any of its components. Modification of the device may result in an unpredictable operation that could cause injury to the patient or operator. Modification of the device voids the warranty.
- The stretcher must only be used in emergency vehicles equipped with a compatible securing system.
- The compatible fastening system in the emergency vehicle must be correctly installed to ensure proper securing of the stretcher. Failure to do so may result in injury to the patient or operator, or damage to the vehicle.
- Practise moving the stretcher between the loading and unloading positions on the rail until you are fully familiar with its operation. Improper use could lead to injury.
- Do not allow untrained assistants to help you operate the rail. Inadequately trained technicians and assistants may cause injury to the patient and to themselves.
- Do not climb on the base of the rail as this may damage the product, unbalance the patient and/or cause injury.
- Avoid placing hands, fingers or feet on the moving parts of both the rail and the stretcher. To avoid injury, when raising and lowering the stretcher, be very careful when placing hands and feet near the moving parts (inlet tray).
- Transporting the stretcher on its side could cause it to tip over and result in damage to the product and injury to the patient. Transporting the stretcher in a lowered position reduces the chances of the stretcher tipping over.
- Always use all the restraining belts to secure the patient on the stretcher. If not properly secured, the patient could fall off the stretcher and suffer an injury.
- Never leave a patient unattended on the stretcher, as they may suffer an injury. When transporting a patient, hold the stretcher firmly in place.
- Side rails are not designed to be used as a means of restraint.
- The TG-880C stretcher is designed to be compatible with the R-450N, R-450N INOX, R-419L and R-800L rails, and with the R-450S INOX, R-419 and R-800G rails for the TG-880S C4 model. It is the responsibility of the operator to ensure that these products work together.
- If the stretcher is occupied by a person of large size, we recommend two operators to handle the stretcher and secure the patient.
- Operators must be able to lift the full weight of the patient, stretcher and additional stretcher elements.
- When cleaning, use appropriate personal protective equipment (protective goggles or respirators) to avoid the risk of inhaling infectious organisms.

- Some cleaning products are corrosive and may damage the product if used incorrectly. If the products described above are used to clean Kartsana equipment, measures must be taken to ensure that the stretchers are rinsed with clean water and thoroughly dried after cleaning. Failure to properly rinse and dry the units may leave a corrosive residue on the surface of the units, which may cause premature corrosion of critical components.

- Failure to properly clean or improperly dispose of a contaminated mattress or other components increases the risk of exposure to blood-borne pathogens and could result in injury to the patient or operator.

PRECAUTIONS

- Changes or modifications to the unit not expressly approved by Kartsana may void the user's authority to operate the system.

- The installation of the stretcher-compatible fixing rail must be carried out by a certified mechanic familiar with the ambulance structure. Consult the vehicle manufacturer before installing the rail and ensure that the installation does not damage or interfere with the vehicle's brake, oxygen or fuel hoses, fuel tank or electrical wiring.

- Before putting the stretcher into operation, remove all obstacles that may interfere and cause injury to the operator or patient.

- When unloading the stretcher from the patient compartment, make sure that the wheels are placed securely on the ground. Failure to do so may result in damage to the product.

- Ensure that the safety belts do not become entangled in the base frame when raising and lowering the stretcher.

- Do not store any objects under the stretcher mattress. Storing objects under the mattress may interfere with the functioning of the stretcher.

- Do not use steam or ultrasonic cleaners to clean the unit.

- The maximum temperature of the water must not exceed 80°C.

- Allow the stretcher to air dry.

- Failure to comply with these instructions may result in some or all of the warranties becoming void.

- A preventive maintenance programme should be established for all Kartsana equipment. Depending on the product's frequency of use, preventive maintenance may need to be carried out more frequently.

- For further information on maintenance, please refer to our technical department.

- Improper maintenance may cause injury or damage to the product. The stretcher should be maintained as described in this manual. Follow the maintenance procedures and use only Kartsana approved parts. The use of unapproved parts or procedures may result in an unpredictable operation and may cause injury and void the warranty.

- Failure to use authorised spare parts, lubricants, etc. may damage the stretcher and void the product warranty.

NOTES

- Loose objects or a build-up of dirt on the patient compartment floor may interfere with the operation of the rail and stretcher securing system. Keep the floor of the patient compartment clean.

- This manual must be considered part of the rail and, as such, must be kept with the product, even if it is sold.

- Kartsana works constantly to improve the design and quality of its products. Therefore, although this manual contains the most current information available, there may be minor discrepancies between the stretcher and this manual. If you have any questions, please contact Kartsana's customer service or technical service area.

- Before installation, Kartsana recommends that a certified technician plans the positioning of the fastening rail inside the emergency vehicle.

2. INSTRUCTION MANUAL

2.1. Technical specifications and dimensions of the TG-880 + MTR stretcher-trolley

Do not damage the areas where the stretcher mechanisms are located, as this may cause malfunction.



Dimensions and measurements in mm.
TG-8804 and MTR-120 models

(*) For the dimensions of each product in the series, please refer to the specifications table in the previous section.


CERTIFIE
10G

kg
MAX WEIGHT
MAX PESO
280 kg


Stretcher with folded legs:


2.1.1. Stretcher-to-trolley locking mechanism

Before locking or unlocking the stretcher, ensure the trolley wheels are braked (see Section 2.3.5. of this manual).


To unlock the stretcher, press and hold the red button as shown in the figure and pull the stretcher from front to rear. Once the intermediate axle is released from the latches, the stretcher can be separated from the trolley.

To lock the stretcher, position it on the rear section of the wings of the TG-880 trolley's longitudinal profile. Once the stretcher is correctly positioned on the trolley, proceed to lock it by sliding the stretcher from rear to front until the red hook engages the stretcher's intermediate axle.

The label indicates that the front wheels are unlocked in the left-hand position and locked in the right-hand position.

2.2. Technical specifications and dimensions of the MTR-120 series stretcher

Do not damage the areas where the stretcher mechanisms are located, as this may cause malfunction.

Dimensions and measurements in mm.
MTR-120 model


(*) For the dimensions of each product in the series, please refer to the specifications table in the previous section.

CERTIFIE
10G

kg
MAX WEIGHT
MAX PESO
280 kg

Note: The front part of the stretcher corresponds to the backrest section of the stretcher.

2.2.1. Operating the stretcher


To minimise the risk of injury to the patient and occupants in the event of an accident, it is advisable to avoid sharp edges and protruding surfaces inside the ambulance, paying special attention to the areas near the stretcher. To avoid potential injury to limbs and other body parts, do not place them along the path of the moving parts of the stretcher.

2.2.2. Removable handles

Extendable handles

To adjust the handle length, press the button on the top and gently pull the handle outwards to the desired length. Continue pulling until the mechanism reaches its stop position, at which point it will be securely locked in place for use.

To return the handle to its original position, press the top button again and gently push the handle inward until it reaches its initial stop. Ensure that the system is properly locked before using the stretcher.


Removable intermediate handles

To remove the two intermediate handles, pull directly on the black knob until it reaches the stop position.

To return the handle to its original position, push the black knob as far as the mechanical stop; a spring system keeps these knobs in this position.

2.2.3. Adjustable backrest section


To adjust the backrest section, operate the lever located directly beneath it and move it to the desired position. Then release the lever to lock the backrest section in the selected position.

2.2.4. Adjustable leg section

By pressing button **A** and manually adjusting the leg section, the leg position can be set to the desired angle. Once the angle is correct and button **A** is released, the clamp is automatically locked in the required position.


Using button **B**, the foot section can be adjusted to the desired angle without altering the leg section. Operate the button in the direction indicated in the figure below and release it once the desired position is reached.

Note: this function is only available on the MRT-120 model

2.2.5. Safety railing

The safety railing has a red handle that locks it in place. By pressing the lever, the rail can be folded down as shown below; when it is manually raised back to its original position, it will lock automatically.

The railing can also be folded sideways by pulling the handles shown below. When the railing is manually moved to its upright position, it will lock automatically.


The safety railings must be in their upright or folded position when loading or unloading the stretcher onto or from its rail.

2.2.6. IV pole

The stretcher is equipped with an IV pole located on the left-hand side. Using knob **A** located at the bottom, secure the IV pole assembly to the side in the desired position.

Loosen knob **B** to raise the IV pole to a vertical position, perpendicular to the stretcher, then tighten the knob again.


The maximum permissible load is 6 kg.


To adjust the IV pole length, press the small positioning button on the inner side of the assembly, then raise the hook to the upper position.

2.2.7. Safety belts

To secure the safety belts on the backrest section, place the buckles (**A**) in the position shown in the figure above. Pass the male buckle (**B**) between the anchors (**A**) until it engages with buckle (**C**). Once secured, centre the assembly and adjust the strap tension at the ends to tighten it.

To secure the patient's abdominal strap, position the female buckle (A) and male buckle (B) as shown in the figure below, then insert the male buckle into the slot of the female buckle. Once secured, centre the assembly and adjust the strap tension to tighten it.

To release the patient, press the button marked 'PRESS'. The securing mechanism will release, allowing the strap buckles to be removed from the locking system.

The patient should be secured with the safety belts at all times while on the stretcher.

2.2.8. Removing and fitting the stretcher belts

Mounting via plate bracket:

Loosen the screws as shown in the figure.

For assembly, follow the reverse procedure: insert and tighten the screws into the nuts located inside the aluminium profile channel.

Position of the belts:

The safety belts located on the backrest section must be positioned over the patient, providing restraint from the shoulders to the abdomen.

The backrest section belts are secured to the profile and must be threaded through the openings provided in the backrest section for correct positioning.

An additional safety belt is provided in the foot section to ensure proper immobilisation.

Securing the safety belts

2.3. Incubator version

2.3.1. Product description and dimensions

The **MT-118IN** stretcher has been specifically designed for the transport of incubators. The product complies with the MDR (EU) Regulation 2017/745. This stretcher has been tested in accordance with the UNE-EN ISO 1865-3 and UNE-EN 1789 standards.


This product has the same measurements and functionalities as the **MTR-120** stretcher, except for the upper part, which does not have a backrest section and leg section, but consists of a metal plate where the incubator is fixed. The main characteristics of the stretcher can be found below.

Dimensions and measurements in mm.
MT-118IN

(*) For the dimensions of each product in the series, please refer to the specifications table in the previous section.

2.4. Technical specifications and dimensions of the TG-880 series trolley

Do not damage the areas where the stretcher mechanisms are located, as this may cause malfunction.

Dimensions and measurements in mm.
TG-880 C2 and TG-880 C4 models



(*) For the dimensions of each product in the series, please refer to the specifications table in the previous section.

CERTIFIE
10G


kg
MAX WEIGHT
MAX PESO
280 kg

Note: The front part of the stretcher corresponds to the backrest section of the stretcher.

Trolley with folded legs:


2.4.1. Rear controls for locking/unlocking the rear wheels

To avoid potential injury to limbs and other body parts, do not place them along the path of the moving parts of the trolley.

Note: The rear part of the trolley corresponds to the foot end of the stretcher.

The black trigger on the rear housing locks or releases the swivel function of the front wheels. See figure.

The label on the housing shows that the front wheels are unlocked in the left position and locked in the right position.

2.4.2. Controls for operating the stretcher legs

Rear controls:

As indicated on the rear housing label, the orange lever operates the folding/unfolding of the front leg, while the blue lever operates the rear leg.

The rear trigger initiates the preliminary folding of the front leg for loading into the ambulance.

To avoid possible injury, do not place any body parts in the path of the moving parts of the trolley.




Front controls:

Engage the intermediate positions of the trolley's front leg (starting from the fully extended position).

Note: To return the leg to the fully extended position, operate the front control and manually lift that section of the trolley so the leg returns automatically. Release the control just before reaching the initial position and check that it locks with an audible click.

To operate the red handle, unfold the movable handle to unlock it, then pull the red control.

Both the red handle and the blue lever located at the front perform the same function.

2.4.3. Intermediate positions of the stretcher trolley

The stretcher has three intermediate height positions relative to the ground, depending on requirements.

To select any of these positions, use the trigger on the backrest section to operate the front leg and the trigger on the leg section to operate the rear leg.


To ensure safe patient handling, the intermediate position must always be operated by two people. Before operating the intermediate positions, secure the patient with the safety belts and raise the side rails.

Hold the stretcher firmly with both hands to prevent dangerous movement.

Intermediate position

1

Operate the front and rear levers simultaneously. Release them at the desired height and check that they are locked.

To return to the previous position, keep both triggers engaged and raise the stretcher using the handles.

Reverse Trendelenburg position

2

Operate the rear lever. Release it at the desired height and check that it is locked.

To return to the previous position, keep the trigger engaged and raise the stretcher using the handles.

Trendelenburg position

3

Operate the front lever. Release it at the desired height and ensure it is locked.

To return to the previous position, keep the trigger engaged and lift the stretcher using the handles.

2.4.4. Brake lever

Operating the levers with the foot locks the rear wheels.

To unlock them, use the tip of your foot to lift the lever.

This special brake prevents both longitudinal movement and swivelling of the trolley.

2.4.5. Mounting the trolley on the ambulance stretcher holder

Move the trolley to the ambulance platform and position the front leading wheels (100 mm diameter) onto the stretcher support. If the platform is fitted with a Kartsana rail, place these leading wheels outside the rails.

Then operate the rear trigger of the front leg to fold it and push the trolley into the stretcher support until it reaches the rear leg. At this point, operate the rear leg control to fold it.

After operating the trigger, you must support part of the trolley's rear weight, so hold it firmly with both hands.

Kartsana specifies that the stretcher trolley must be secured to the ambulance using compatible Kartsana rails (see table), specifically designed for this purpose.

MODEL	TG-8802	TG-8804
Compatible securing systems	R-419L R-450N R-450N INOX R-800L	

Once the described process has been completed, the trolley can be fully inserted into the rail until it is locked by the front and rear anchors. The Kartsana rail secures the stretcher trolley by means of two front anchors and one rear anchor. These anchors prevent the stretcher trolley from tipping over in the event of an accident.

If the specified rails are not used, Kartsana does not guarantee the functionality or safety of a system that has not been tested by the company.

3. CARE AND MAINTENANCE OF THE PRODUCT

3.1 Contact details

If you have any questions regarding the use, maintenance or installation of a KARTSANA product, please contact our Technical Service Area (TSA).

Telephone: +34 93 715 86 72
Email: sat@kartsana.com
Address: Narcís Monturiol, 34
08192 Sant Quirze del Vallès,
Barcelona (Spain)

Any serious incident related to the product must be reported to the manufacturer and to the competent authority of the Member State where the user and/or patient are located.

3.2 Environmental conditions

Both for use and storage:

Recommended temperature: from -10°C to +60°C

Relative humidity: 5–95%

3.3 Cleaning

To clean the stretcher, apply suitable cleaning products that do not damage the material surfaces using a cloth, rag or similar. If the surface is lubricated, re-lubricate it after cleaning.

Please use cleaning agents according to the manufacturer's instructions.

The manufacturer will not be liable for any damage or anomaly caused by a cleaning product that may damage the surface of the stretcher materials.

For reasons of hygiene and preservation of the product components, the manufacturer recommends cleaning after each use.

3.4. Maintenance

3.4.1 Precautionary maintenance

Correct and regular maintenance guarantees the durability of the product.

Creating a maintenance plan, including regular reviews and establishing a responsible person to carry them out, is recommended.

The person in charge must adhere to the following **requirements**:

- Have technical knowledge related to the product and the maintenance procedures described in these instructions.
- Have the support of qualified, trained and knowledgeable staff who are familiar with the operation of the product.
- Use original or manufacturer-approved components, spare parts and accessories.
- In accordance with the instructions of the European Council Directive 93/42/EEC which establishes the obligation of the purchaser to provide on request the above mentioned after-sales care record for the purpose of product traceability, records of all product maintenance operations will be kept.

Check for proper operation **before use**. Prior to service, the following should be checked:

- Functionality.
- Fastenings and hardware.
- Parts that move, such as wheels and belts.
- Sensitivity of springs.

If the device does not appear to be suitable for correct and safe use, **it should be taken out of service** until the device is repaired or restored to full working order.

Do not modify the structure of the device for any reason, as this may cause serious injury to patients and/or operators.

Lubrication intervals of moving parts:

- 1–30 services per month: every 3 months.
- 30–50 services per month: every 2 months.
- More than 150 services per month: every month.

Check and lubricate the inside of the leg channel every two months.

All maintenance other than lubrication, tightening of nuts and bolts and ordinary cleaning must be carried out by KARTSANA or an authorised service centre.

3.4.2. Service maintenance

The person responsible for inspecting the product must meet the following basic requirements:

- ◆ Sufficient knowledge of the product, its technical/constructional characteristics, control points and final tests, packaging, preservation and handling.
- ◆ Knowledge of all product functions and any possible risks, malfunctions or failures of the product.
- ◆ Possession of all necessary tools to carry out any technical service or minor repairs.
- ◆ Possession (or ability to procure) original manufacturer's or manufacturer's authorised spare parts.
- ◆ Use or support of specialised technical personnel trained by the manufacturer to service the product in question.
- ◆ Keeping a record of all maintenance operations carried out on the product, in accordance with the instructions of the European Council Directive 93/42/EEC which establishes the obligation of the purchaser to provide on request the above mentioned after-sales care record for the purpose of product traceability.

3.5 Transport and storage

Before transporting the device, make sure that it is properly packaged and secured to prevent damage during transport.



Keep the original packaging for possible further transport. Any damage caused to the device during transport is not covered by the warranty. Repairs and replacements of damaged parts are the responsibility of the customer. The device should be stored in a dry place.

When storing, do not place any heavy objects on top of the product. It must not be used as a support for any kind of object.

3.6 Disposal

When the product is no longer in a usable condition, do not throw it away but develop a take-back strategy for recycling by an authorised dealer. The packaging should be treated as reusable waste. The metals are fully recyclable. The plastic products should be treated as recycling material. Waste management should be as established in each country. If you have any questions, please ask your municipal administration.

4. SPARE PARTS

To consult the spare parts for this product, you can visit the **following link**:

https://www.kartsana.com/global/Spare_parts_MTR120_TG880.pdf

If you have any questions, please contact our technical service.

5. STRETCHER ACCESSORIES

C-003 Oxygen bottle holder

C-011 Jupiter foldable bariatric kit

C-041 Corpuls Jupiter/241 support

C-063 Leg folding system


C-079
C-084 Bariatric click for Jupiter

C-102 Instrument table 880-241

C-112 Superbrava rail kit for mechanical stretchers

C-114 Standard cushion


C-116 SCH monitor mount for mechanical stretchers

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