
# Serie

# 241

## TGX-241 (M/L/XL)

**KARTSANA**

Life-moving innovation

C/ Narcís Monturiol, 34  
08192 Sant Quirze del Vallès  
Barcelona (España)  
+34 93 715 86 72  
info@kartsana.com  
www.kartsana.com

Manual de instrucciones

**Camilla mecánica TGX-241**

Instruction manual

**Mechanical stretcher TGX-241**

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**Brancard mécanique TGX-241**

ESPAÑOL

ENGLISH

FRANÇAIS

# KARTSANA

Life-moving innovation

C/ Narcís Monturiol, 34  
08192 Sant Quirze del Vallès  
BARCELONA  
TEL. +34 93 715-49-24  
[info@kartsana.com](mailto:info@kartsana.com)  
[www.kartsana.com](http://www.kartsana.com)

## Serie 241

# CAMILLA MECÁNICA

Models TGX-241 M / L / XL

MANUAL DE INSTRUCCIONES

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Models TGX-241 M / L / XL

INSTRUCTIONS MANUAL

Please read the instruction manual before use and keep for future reference.

CERTIFIED  
**10G**  
EN-1789

CE

## INDEX

---

<b>1.- INTRODUCTION</b>	<b>3</b>
1.1.- FOREWORD:	3
1.2.- RESPONSIBILITY AND WARRANTY	3
1.3.- SPECIFICATIONS	3
1.4.- ATTENTION	4
1.5.- SUMMARY OF SAFETY PRECAUTIONS	4
<b>2.- INSTRUCTION MANUAL</b>	<b>6</b>
2.1.- TECHNICAL CHARACTERISTICS	6
2.2.- OPERATION AND HANDLING OF THE STRETCHER	7
2.2.1.- EXTENSIBLE HANDLES	7
2.2.2.- ADJUSTABLE RECLINING BACKREST	7
2.2.3.- ADJUSTABLE LEG REST	9
2.2.4.- SAFETY RAIL (ALL MODELS)	9
2.2.5.- IV DRIP STAND	11
2.2.6.- HOW TO SECURE AND FASTEN THE SAFETY BELTS	12
2.2.7.- FASTENING THE MIDDLE SAFETY BELT	12
2.2.8.- REMOVAL AND RE-ATTACHMENT OF SAFETY BELTS TO THE STRETCHER	13
2.2.9.- LEG POSITION CONTROLS	15
2.2.10.- INTERMEDIATE STRETCHER POSITION, TRENDELENBURG & COUNTER-TRENDELENBURG	16
2.2.11.- BRAKE LEVER	18
2.2.12.- MOUNTING THE STRETCHER ON ITS RAIL FIXATION SYSTEM IN THE AMBULANCE	19
<b>3.- CLEANING &amp; DISINFECTION</b>	<b>20</b>
<b>4.- MAINTENANCE</b>	<b>20</b>
<b>5.- EXPLODED VIEWS</b>	<b>21</b>

# OPERATION AND HANDLING OF TGX-241 M, L & XL STRETCHERS

## 1.- INTRODUCTION

### 1.1.- Foreword:

The TGX-241 stretchers are specially designed for the rescue and transportation of patients.

The TGX-241 stretchers are and have been tested in accordance with the UNE-EN 1865 and UNE-EN 1789 standards.

All the information regarding treatment, disinfection and maintenance is indicated, taking into consideration our current experience and know-how.

Certain technical modifications to the stretcher have been placed under reserve for the purpose of making improvements to the product.

### 1.2.- Responsibility and warranty

The stretcher must be checked at the time of its delivery to the assistance organisation. All its functions must be explained in great detail. The assistance organisation must take charge of teaching all its employees how to use it correctly.

The warranty is valid for two years following the date indicated on the purchase invoice.

The warranty will not cover the fault if that fault is due to incorrect installation, mistreatment or improper use of the stretcher. All repairs must be made by a technical service authorised by KARTSANA or its respective representative.

### 1.3.- Specifications

MODELS	TGX-241 M, L, XL
Maximum working load	250 Kg*
Standards	UNE-EN-1865-1 EN-1789+A1
Total length (non extended, extended)	1953 / 2254 mm
Total width	590 mm
Loading Height	595 mm (M) 642 mm (L) 686mm (XL)
Weight	42,6 Kg (M) 43,8 Kg (L) 45 Kg (XL)
Recommended operatives required to load / unload an occupied stretcher	1 operative, or 2 if the stretcher were occupied by a larger sized patient.
Compatible fixation rail systems	R-419 / R-450-S Inox / R-800 & R-900

\*Stretcher not valid for bariatric or paediatric transport.

#### 1.4.- Attention

This symbol indicates important safety measures for the correct use of the stretcher, in order to prevent accidents. The warnings alert the reader about a situation which, if not avoided, could result in death or serious injury. The cautions alert the reader of a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient, or damage to the equipment or other property. This includes special care necessary for the safe and effective use of the device and to avoid damage that may occur as a result of use or misuse.

#### 1.5.- Summary of safety precautions

Please read the warnings and cautions listed on these pages carefully. Repair and maintenance services must be carried out by qualified personnel only.

#### WARNINGS:

- Improper use of the stretcher can cause injury to the patient or operator. Operate the stretcher only as described in this manual.
- Do not modify the stretcher or any components of the stretcher. Modifying the product can cause unpredictable operation resulting in injury to the patient or operator. Modifying the product also voids its warranty.
- Any emergency vehicle to be used with this stretcher must have the compatible fixation system installed.
- Have the compatible fixation rail installed by a certified mechanic. Improper rail installation can cause injury to the patient or operator and/or damage to the stretcher. Verify that the stretcher legs lock into the load position before without contact with the fender of the vehicle. Failure to properly lock the stretcher height into position can cause injury to the patient or operator and/or damage to the stretcher
- Practice changing height positions and loading the stretcher until operation of the product is fully understood. Improper use can cause injury.
- Do not allow untrained assistants to assist in the operation of the stretcher. Untrained technicians/assistants can cause injury to the patient or themselves.
- Grasping the stretcher improperly can cause injury. Keep hands, fingers and feet away from moving parts. To avoid injury, use extreme caution when placing your hands and feet near the base tubes while raising and lowering the stretcher.
- Always use all restraint straps to secure the patient on the stretcher. An unrestrained patient may fall from the stretcher and be injured.
- Never leave a patient unattended on the stretcher or injury could result. Hold the stretcher securely while a patient is on the product.
- Side rails are not intended to serve as a patient restraint device.
- High obstacles such as curbs, steps or rough terrain can cause the stretcher to tip, possibly causing injury to the patient or operator.
- The TGX-241 stretchers are compatible with the following KARTSANA fixation rail systems: R-419 / R-450-S Inox / R-800 & R-900; it is the operator's responsibility that these products work together.
- Two operators must be present when the stretcher is occupied.
- Operators must be able to lift the total weight of the patient, stretcher and any items on the stretcher.
- Never install or use a wheel lock on a stretcher with excessively worn wheels. Installing or using a wheel lock on a wheel with a diameter inferior to 200 mm could compromise the effectiveness of the wheel lock, possibly resulting in injury to the patient or operator and/or damage to the stretcher or other equipment.
- When cleaning, use any appropriate personal safety equipment (goggles, respirator, etc.) to avoid the risk of inhaling harmful or infectious substances.

# KARTSANA

Life-moving innovation

## OPERATION AND HANDLING OF TGX-241 M, L & XL STRETCHERS

- Some cleaning products are corrosive in nature and may cause damage to the product if used improperly. If the products described above are used to clean KARTSANA equipment, measures must be taken to ensure the stretchers are wiped with clean water and thoroughly dried following cleaning. Failure to properly rinse and dry the stretchers will leave a corrosive residue on the surface of the stretchers, possibly causing premature corrosion of critical components.
- Failure to properly clean or dispose of a contaminated mattress or other stretcher components will increase the propagation risk of bloodborne pathogens and may cause injury to the patient or operator.

### CAUTIONS:

- Changes or modifications to the unit not expressly approved by KARTSANA could void the user's authority to operate the system.
- Installation of the fixation rail system should be done by a certified mechanic familiar with ambulance vehicle construction. Consult the vehicle manufacturer before installing the rail system and be sure that the installation does not damage or interfere with the brake lines, oxygen lines, fuel lines, the fuel tank or electrical wiring of the vehicle.
- Before operating the stretcher, clear any obstacles that may interfere and cause injury to the operator or patient.
- When unloading the stretcher from the patient compartment, ensure that the caster wheels are safely set on the ground or damage to the product may occur.
- Wheel lock(s) are only intended to help prevent the stretcher from rolling away while unattended, and to aid in patient transfer. A wheel lock may not provide sufficient resistance on all surfaces or under loads.
- Ensure that the restraints are not entangled in the base frame when raising and lowering the stretcher.
- Do not store items under the stretcher mattress. Storing items under the mattress can interfere with the operation of the stretcher.
- Do not steam clean or ultrasonically clean the unit.
- Maximum water temperature should not exceed 80°C.
- Allow stretcher to air dry.
- Towel dry all casters and interface points.
- Failure to comply with these instructions may invalidate any or all warranties.
- Improper maintenance can cause injury or damage to the product. Maintain the stretcher as described in this manual. Use only KARTSANA approved parts and maintenance procedures. Using unapproved parts and procedures could cause unpredictable functioning and/or injury and will void the product warranty.
- Failure to use authorized parts, lubricants, etc. could cause damage to the stretcher and will void the warranty of the product.

### NOTES:

- Loose items or debris on the patient compartment floor can interfere with the operation of the safety hook and stretcher fastener. Keep the patient compartment floor clear
- This manual should be considered a permanent part of the stretcher and should remain with the product even if the stretcher is subsequently sold.
- KARTSANA continually seeks advancements in product design and quality. Therefore, while this manual contains the most current product information available at the time of printing, there may be minor discrepancies between your stretcher and this manual
- KARTSANA recommends that, prior to installation, the certified mechanic should plan the placement of the fixation system in the rear of the vehicle.

## 2.- INSTRUCTION MANUAL

The stretcher is adapted for mounting on the Kartsana R-419, R-450-S, R-700 and R-800 rails, based on the UNE-EN 1789 standard.

To reduce the risk of injuries to the patient and accompanying persons in the event of an accident, it is advisable to avoid sharp edges and projecting parts inside the ambulance, paying special attention to the areas nearest to the stretcher.


Make sure that no damage is caused to the areas containing the stretcher mechanisms, to prevent them from malfunctioning.

### 2.1.- Technical characteristics

Aproximate stretcher weight: 45 Kg.

Maximum load capacity: 250 Kg.

Maximum measurements shown below:

### 2.2.- Operation and handling of the stretcher


To prevent injuries to limbs and other body parts, make sure they are not within the range of the moving parts of the stretcher.

Note: The “front” part of the stretcher refers to the backrest area where the head and torso of the patient would rest.

#### 2.2.1.- Extensible handles


To extend the handles, press down on the small red button and pull them out as far as they go.

To return the handles to their original position, press the button and slide them in again.

#### 2.2.2.- Adjustable reclining backrest

To adjust backrest inclination, lift the lever behind the patient's head area and move the backrest to the desired position, then release the lever to lock the position.

**- Foldable front**


# KARTSANA

Life-moving innovation

## OPERATION AND HANDLING OF TGX-241 M, L & XL STRETCHERS

As shown in the following image, the entire front of the stretcher can be folded to a vertical position to allow the stretcher to manoeuvre in tight spots:



This is done by pressing the lever circled in red and lifting the handles to the folded position.

### 2.2.3.- Adjustable leg rest

#### Button for adjusting leg rest inclination


The required inclination of the leg rest can be obtained by pressing button A and moving the leg rest manually. Once the correct inclination has been reached, release button A to automatically lock the leg rest.


The adhesive to the right indicates the position of the button and its function

### 2.2.4.- Safety rail (all models)

The safety rail has a red bar that locks the rail into place. To lower the rail, just press the bar as shown in the figure below. To lock it in the safety position again, move it to that position manually, and it will automatically lock into place.

The safety rail can also be folded outwards. To fold the rail, just press the 2 red levers as shown in the figure below. To lock it in the safety position again, move it to that position manually, and it will automatically lock into place.


**ATTENTION: The safety rails must be in their vertical or lowered position to load or unload the stretcher onto its rail.**


### 2.2.5.- IV drip stand

The stretcher's IV stand rests along the left side of the "back" of the stretcher.

By loosening handle A the IV drip stand can be removed to reposition it on either side of the stretcher; loosening handle B, the IV drip stand can be moved to its upright position when needed, locking it in place by tightening handle B.


The length of the IV drip stand can be adjusted by pressing a tiny metal button along the inside of the shaft, and pulling up on the hook till the shaft extends as far as it will go.

The maximum load the IV drip stand can hold is 6 Kg.

The maximum load is indicated on a sticker like this one:

PESO MÁXIMO:  
**6 Kg**  
 MAXIMUM LOAD:

### 2.2.6.- How to secure and fasten the safety belts


Place the clasps (C) in the position shown in the figure above. Pass the male part of the buckle (D) between the openings (C) until it is inserted into the female part B. Once the assembly is anchored in place, centre and tauten the assembly, adjusting the tension of the belt at ends A and/or E.

NOTE: The foot belt is not fitted with the C openings.

### 2.2.7.- Fastening the middle safety belt

It is advisable for this belt to be attached to the stretcher to ensure the complete safety of the patients and guarantee compliance with the UNE-EN 1865 standard test that is performed.


Put the female part of the buckle A and male part B in the position shown in the figure below and insert the latter into the slot in the former. Once the assembly is firmly secured, centre and tauten it, adjusting the tension of the belt at end C.

It is advisable to ensure that while the patient is secured by the safety belts at all times while on the stretcher.

### 2.2.8.- Removal and re-attachment of safety belts to the stretcher

The safety belts are attached to the stretcher as so:

**A:** The straps of the chest-belts are threaded through the slot and tied around the crosspiece using the opening around it, as shown above.

**B:** The middle belts are tied to the side bar using the available slot.


**C:** The leg belt is threaded through the hole along the side of the leg rest and tied around the tube as shown.

### 2.2.8.1.- Safety Belt Positions

**2.2.8.1.1.-** Middle and leg belts are positioned at the approximate distances shown in the figure below:



**2.2.8.1.2.-** The chest belts are threaded through the slots in the backrest and tied to the crosspiece as shown on the previous page.



## 2.2.9.- Leg position controls

### 2.2.9.1.- Back controls

- Blue handle (right): Allows folding / unfolding of the stretcher's front legs.


- Orange handle (left): Allows folding / unfolding and intermediate positioning of the stretcher's back legs.

- Note: To collapse legs into their folded position, release the handle just before they reach their folded position. To set the legs in their intermediate positions, release the handles as soon as the legs begin to fold; if the position the legs lock into is not the desired one, repeat the process until the required position is achieved. A “click” sound can be heard each time the legs lock into a position. Listen for this sound to be sure they are correctly locked.

### 2.2.9.2.- Lever for front leg intermediate positioning

- 
- The red lever under the left stretcher handle accesses the intermediate positions of the front legs, starting from the unfolded (upright) position.
- Note: To return the front legs to their original unfolded position, activate the red lever and then raise the front end of the stretcher manually so that the leg falls automatically into place. Just before the leg reaches upright position, release the lever so that when the leg reaches that position, the click can be heard, indicating the leg is securely locked in place.


### 2.2.10.- Intermediate stretcher position, Trendelenburg & Counter-Trendelenburg

- The stretcher has three intermediate positions at different heights from the ground.
- To ensure the safe treatment of the patient, the intermediate position should always be activated by two people. In addition, the safety belts must be placed on the patient and the rails raised before activating the intermediate positions.
- To avoid dangerous unpredictable movements of the stretcher, be sure to hold on to the handles firmly with both hands.
- Either end of the stretcher can be brought to intermediate height using the aforementioned controls.
- 1st position: Intermediate position


- It is achieved by activating the intermediate position lever (red) located on the right-hand side of the stretcher. This releases the front legs down to a mechanical stop.
- Then, by pressing and immediately releasing the orange handle, the back legs will fall into first intermediate position.
- To return the stretcher to upright position, slightly lift the stretcher and using one foot gently help the front legs to rise.

- 2nd position: Trendelenburg position


- It is achieved by activating the intermediate position lever (red) located on the right-hand side of the stretcher. This releases the front legs up to a mechanical stop.

- 3d position: counter-

Trendelenburg position

- It is positioned by activating the rear control (orange) once until it reaches the first position..


### 2.2.11.- Brake lever

- Activate the red lever with your foot to lock the rear wheels. This blocks both the longitudinal and the rotating movement of the wheel.

- - How to lock / unlock swivel casters
- The green bolt is used to lock or unlock wheel rotation, as shown in the following figure. When in intermediate position, the green bolt must be in locked position at all times, as the sticker indicates:

### 2.2.12.- Mounting the stretcher on its rail fixation system in the ambulance


- KARTSANA recommends using the KARTSANA rails to secure the stretcher to the ambulance, since these have been specially designed for that purpose.

- Guide the stretcher to the ambulance platform and mount its front wheels (the smaller 10mm diameter wheels), into the the stretcher dock; if the platform has a KARTSANA rail, place the front wheels just before the rail. Activate the blue handle to fold the front legs and slide the stretcher into the stretcher dock up to the rear legs. At this point, activate the Orange handle to allow the rear legs to fold.

- After activating the orange handle, part of the weight of the rear end of the stretcher must be supported, and so it will be necessary to hold it firmly with both hands.

- Once the above process has been executed in full, the stretcher can be inserted completely into the stretcher dock, until it is locked in place at the front and rear. The KARTSANA rail fixation system anchors the stretcher at the front by means of two anchoring points at the front and one anchoring point at the rear (see the figures below). These anchoring points prevent the stretcher from overturning in the event of an accident. In the event of not using the KARTSANA rails, it is advisable that the front part of the stretcher cart be anchored at the end of the bar marked 12 in the exploded diagram of the cart that is included in this manual.


### 3.- CLEANING & DISINFECTION


- When disinfecting the stretcher, use products that will not damage the surface materials and wipe it down with a cloth or similar. If the surface to clean is greased, re-grease it after disinfecting it.
- All disinfectants must be used in accordance with the manufacturer's instructions.
- The manufacturer will not be held responsible for any anomaly or damage caused by using a cleaning product that could damage the surface of the rail materials.
- To ensure the hygiene and proper conservation of product components, the manufacturer recommends that they be disinfected after each use.

### 4.- MAINTENANCE

- Check and grease the interior of the guide (shown below) every six months.
- The stretcher wheels and the locking systems must also be greased. We recommend that maintenance be performed by a specialist from our company every six months. During this inspection, a general check will be made to replace any parts that could be damaged or worn due to use.

## 5.- SPARE PARTS

For spare parts for this product, please visit the following link:

[https://www.kartsana.com/global/Spare\\_parts\\_TGX241.pdf](https://www.kartsana.com/global/Spare_parts_TGX241.pdf)

If you have any questions or queries, please contact our technical service.

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