



DATA SHEET

MT-118 / MT-120 / MTR-120

JUPITER TROLLEY STRETCHER

MECHANICAL



Aluminium



Guardrails
*MT-120,
MTR-120



Articulated
*MTR-120

CERTIFIED
10G



MT-118:

- ▷ Adjustable headrest.
- ▷ Safety belt kit for abdomen and legs.

MT-120:


- ▷ Adjustable headrest.
- ▷ Laterally and longitudinally collapsible guardrails.
- ▷ Safety belt kit for shoulders, abdomen and legs.
- ▷ Aluminium tray for performing cardiac massages.

MTR-120:

- ▷ Adjustable headrest.
- ▷ Laterally and longitudinally collapsible guardrails.
- ▷ Safety belt kit for shoulders, abdomen and legs.
- ▷ Aluminium tray for performing cardiac massages.
- ▷ Double-jointed leg rest.

DIMENSIONS AND WEIGHT

	MT-118	MT-120	MTR-120
H Height (mm)	150	150	150
W Width (mm)	570	570	570
L Minimum/maximum length (mm)	1962	1962	1962
Weight (kg)	16,5	17,6	20,9



RELATED PRODUCTS



TG-8802



TG-8804



TG-880 C4 IN



KARTSANA

Life-moving innovation

Narcís Monturiol, 34
08192 · Sant Quirze del Vallès
Barcelona (Spain)

T.: +34 93 715 86 72
info@kartsana.com
www.kartsana.com