



DATA SHEET

# R-419 / R-419E / R-419L

RAILS

MECHANICAL



Aluminium and steel



Ramp



Backboard stretcher drawer


CERTIFIED 10G



- ▷ **R-419:** compatible with 241 Series and JUPITER Series TG-880SC4 model
- ▷ **R-419 E:** compatible with Series 241 and JUPITER Series TG-880SC4 model
- ▷ **R-419 L:** compatible with JUPITER Series except TG-880SC4 model
- ▷ Rail made mainly of aluminum and steel
- ▷ 3 solid anchor points for the stretcher
- ▷ Manual lock / unlock with the input tray
- ▷ Can accommodate spinal or paddle stretcher
- ▷ Fixing to the ambulance with 6 screws DIN 912 M8 -10.9 and 2 screws DIN 7991 M6 - 10.9
- ▷ 10G certification.

## PRODUCT SPECIFICATIONS

	R-419	R-419E	R-419L
<b>H</b> Height (mm)	377	394	347
<b>L</b> Length closed (mm)	2279	2296	2425
<b>L2</b> Length open (mm)	1902	1902	2078
<b>W</b> Width (mm)	623	513	623
Weight (kg)	46	45,5	50



## RELATED PRODUCTS



TG-241



TGE-241



TG-880 C4



KARTSANA

Life-moving innovation

Narcís Monturiol, 34  
08192 - Sant Quirze del Vallès  
Barcelona (Spain)

T.: +34 93 715 86 72  
info@kartsana.com  
www.kartsana.com