



DATA SHEET

R-800G / R-800L / R-800E

RAILS

MECHANICAL



Aluminium and steel



Ramp


CERTIFIED
10G



- ▷ **R-800 G:** compatible with Series 241 and JUPITER Series TG-880SC4 model.
- ▷ **R-800 L:** compatible with JUPITER Series stretchers, except TG-880SC4 model.
- ▷ Rail made mainly of aluminum and stainless steel.
- ▷ 3 solid anchor points for the stretcher.
- ▷ Manual lock / unlock with the input tray.
- ▷ Fixing to the ambulance with 6 screws DIN 912 M8 -10.9 and 1 screw DIN 7991 M8 - 10.9.
- ▷ 10G certification.

PRODUCT SPECIFICATIONS

	R-800G	R-800L	R-800E
H Height (mm)	112	112	108,5
L Length closed (mm)	1869	1959	1857
L2 Length open (mm)	2179	2258	2264
W Width (mm)	623	623	623
Weight (kg)	28,5	42	28,5



RELATED PRODUCTS



TG-241



TGE-241



TG-880 C2



KARTSANA

Life-moving innovation

Narcís Monturiol, 34
08192 · Sant Quirze del Vallès
Barcelona (Spain)

T.: +34 93 715 86 72
info@kartsana.com
www.kartsana.com