



DATA SHEET
TG-1000UK

CHEROKEE SERIES

ELECTRIC EQUIPMENT



Aluminium and stainless steel



Electric



CERTIFIED
10G






- ▷ Adjustable headrest with shock absorber.
- ▷ Double-jointed leg rest with 2 shock absorbers.
- ▷ 4 Ø160 mm rotating wheels (2 lockable wheels + 2 with brakes).
- ▷ Guardrails that collapse longitudinally and fold outwards.
- ▷ Folding IV pole.
- ▷ Three-point safety belt kit.
- ▷ Manual operation system.

PRODUCT SPECIFICATIONS

TG-1000UK

H	Max. height (mm)	948
HC	Max. adjustable height (mm)	900
H2	Min. height (mm)	468
L	Length (mm)	1970
L2	Length closed (mm)	1607
L3	Length open (mm)	2120
W	Width (mm)	567
	Load capacity (kg)	300
	Weight (kg)	69



KARTSANA

Life-moving innovation

Narcís Monturiol, 34
08192 · Sant Quirze del Vallès
Barcelona (Spain)

T.: +34 93 715 86 72
info@kartsana.com
www.kartsana.com