



DATA SHEET
TGM-241

241 COMPACT SERIES

MECHANICAL



Aluminium



Reclinable



250 Kg



Laterally and longitudinally collapsible guardrails



Rotating wheels



Leg rest with shock absorber


CERTIFIED
10G



- ▷ Shock-assisted adjustable head
- ▷ Shock-assisted single-joint leg
- ▷ Trendelenburg, anti-Trendelenburg and intermediate position
- ▷ 4 Ø200 swivel wheels (rear lockable)
- ▷ 4 heights available to adapt to the ramp of the rail of the ambulance (S, M, L and XL)
- ▷ Side and longitudinal folding rails
- ▷ Foldable IV pole
- ▷ Belt kit for shoulders, abdomen and legs (3 points)
- ▷ Aluminum tray for cardiac massage
- ▷ 4 telescopic grips
- ▷ Fireproof welded mattress

DIMENSIONS AND WEIGHT

	S	M	L	XL
H Height (mm)	800	850	900	952
HC Height from the floor to the leading wheel (mm)	560	595	660	681
W Width (mm)	571	571	571	571
L Length (mm)	1951	1951	1951	1951
Weight (kg)	43	43,4	43,8	44,2



RELATED PRODUCTS



R-419



R-800E



KARTSANA

Life-moving innovation

Narcís Monturiol, 34
08192 · Sant Quirze del Vallès
Barcelona (Spain)

T.: +34 93 715 86 72
info@kartsana.com
www.kartsana.com