

KARTSANA


Life-moving innovation



241 COMPACT SERIES



241 COMPACT SERIES



Ability to adapt

A series featuring four models that stand out for their high level of versatility in terms of movement and height.



01

Laterally and longitudinally collapsible guardrails

Easy to transfer patients.



02

Load capacity of 250 kg



03

Double the wheels, double the ease of loading

Two leading wheels (only on the TG model).



04

Safe and easy

Front wheel locking and unlocking (only on the TG model).



05

Wheel lock



One series, five models


		TG	TGM	TGX	TGEM	TGEXL
Stretcher mainly made of aluminium		✓	✓	✓	✓	✓
Adjustable headrest with shock absorber		✓	✓	✓	✓	✓
Leg rest	Single-jointed with shock absorber	✓	✓	✓		
	Fixed				✓	✓
Front wheels	Lockable rotating	✓				
	Rotating		✓		✓	
	Fixed			✓		✓
Rear wheels	Lockable rotating	✓	✓	✓	✓	✓
Leading wheel	Twin	✓				
	Single		✓	✓	✓	✓
3 positions		✓	✓	✓	✓	✓
4 telescopic handles		✓	✓	✓	✓	✓
Laterally and longitudinally collapsible guardrails		✓	✓	✓	✓	✓
IV post		✓	✓	✓	✓	✓
Safety belt kit for shoulders, abdomen and legs (3-point)		✓	✓	✓	✓	✓
Aluminium tray for performing cardiac massages		✓	✓	✓	✓	✓
Fire-resistant welded mattress		✓	✓	✓	✓	✓

Trendelenburg and reverse Trendelenburg

All models in the 241 series fare reclinable.



Adaptable to the height of the ambulance rail ramp



* Recommended loading height, at which the leading wheel makes contact halfway along the rail ramp.

4 available heights



Distance from the floor to the leading wheel to adapt to the rail ramp on the ambulance.

Product certification

All KARTSANA products comply with the EN 1865 and 1789 standards.



Compatibility with rails


Rails	241 COMPACT SERIES
R-419	✓
R-450S STAINLESS STEEL	✓
R-800G	✓

TG-241

DIMENSIONS AND WEIGHT

	S	M	L	XL
H Height (mm)	810	860	904	
HC Height from floor to leading wheel (mm)	530	585	639	-
W Width (mm)	590	590	590	-
L Length (mm)	1951	1951	1951	-
Weight (kg) *	50	51.7	54	-

* Weight with mattress and belts.



Aluminium



Reclinable



250 kg



Laterally and longitudinally collapsible guardrails



Leg rest with shock absorber



Rotating wheels



Lockable wheels



Twin wheels on leading edge

CERTIFIED
10G



TGX-241

DIMENSIONS AND WEIGHT

	S	M	L	XL
H Height (mm)	-	860	909	953
HC Height from floor to leading wheel (mm)	-	590	642	686
W Width (mm)	-	590	590	590
L Length (mm)	-	1953	1953	1953
L2 Length folded (mm)	-	1630	1630	1571
Weight (kg) *	-	49.5	51.5	53

* Weight with mattress and belts.



Aluminium



Reclinable



250 kg




Longitudinally collapsible guardrails



Leg rest with shock absorber




Foldable



TGM-241

DIMENSIONS AND WEIGHT


	S	M	L	XL
H Height (mm)	800	850	910	950
HC Height from the floor to the leading wheel (mm)	570	595	650	705
W Width (mm)	595	595	590	590
L Length (mm)	1951	1951	1951	1951
Weight (kg)	47	47.5	50.6	51.1



TGE-241

DIMENSIONS AND WEIGHT

	S	M	L	XL
H Height (mm)	-	860	-	950
HC Height from the floor to the leading wheel (mm)	-	585	-	700
W Width (mm)	-	590	-	590
L Length (mm)	-	1951	-	1951
Weight (kg)	-	46.5	-	44.1



TGE-241M





KARTSANA S.L.

Narcís Monturiol, 34
08192 · Sant Quirze del Vallès
Barcelona (Spain)
T.: +34 93 715 86 72
E.: info@kartsana.com
www.kartsana.com

07-2022

© 2022 · Kartsana. All rights reserved.