

KARTSANA

Life-moving innovation



INCUBATOR STRETCHERS



INCUBATOR STRETCHERS

INCUBATOR

Safety and care while transporting the most delicate

The widest range for paediatrics offers smooth movement and high levels of autonomy.



TG-1000 IN

- ▷ Stretcher compatible with incubators for paediatric services.
- ▷ Adjustable working height.
- ▷ 4 ø160 mm rotating wheels (2 locking wheels + 2 with brakes).
- ▷ Hydraulic system unlocking button in case of fault.
- ▷ Safety system for avoiding collapse.



Aluminium



Electric



300 kg unassisted



CERTIFIED 10G




TUV



DIMENSIONS AND WEIGHT

		TG-1000-IN
H	Max. height (mm)	1300
H2	Min. height (mm)	350
L	Length (mm)	1972
W	Width (mm)	574
Weight (kg)		61



* Upper C-004 tray not included in reference TG-1000 IN.

01 300 kg unassisted

Stretcher with a hydraulic elevation system which can lift loads of up to 300 kg unassisted.



02 Complete compatibility

Base for installing any model incubator.



03 Unlocking pedal

Wheels are unlocked using a single pedal to give the operator maximum control.



04 Control panel

New and improved control panel, now more intuitive, autonomous and functional.



05 Emergency lights

For increased safety.



06 Maximum autonomy

Features 2 batteries to guaranteed operation.



TG-880 C4 IN



Aluminium



4 lockable wheels



280 kg




CERTIFIED
10G

DIMENSIONS AND WEIGHT

TG-880 C4 IN

H Height (mm)	986
HC Height from floor to leading wheel (mm)	720
W Width (mm)	595
L Length (mm)	1974
Weight (kg)	55



Compatibility with rails

Rails	TG-880 C4 IN
R-419L	✓
R-800L	✓
R-450N	✓

MT-118 IN

- ▷ Compatible with the JUPITER Trolley Stretcher Series.
- ▷ Can be adapted to all models of incubators.
- ▷ 4 telescopic handles + 2 side handles.



Aluminium



CERTIFIED
10G



DIMENSIONS AND WEIGHT

MT-118 IN

H Height (mm)	150
W Width (mm)	590
L Length (mm)	1962
Weight (kg)	32







KARTSANA S.L.

Narcís Monturiol, 34
08192 · Sant Quirze del Vallès
Barcelona (Spain)
T.: +34 93 715 86 72
E.: info@kartsana.com
www.kartsana.com

07-2022

© 2022 · Kartsana. All rights reserved.