

KARTSANA

Life-moving innovation



JUPITER TROLLEY STRETCHER



JUPITER

TROLLEY STRETCHER

JUPITER

Double the functions, double the comfort

Very easy to assemble, makes transporting patients easier by adapting to different spaces.





01

6 telescopic handles

Also features 2 extra handles on the sides (MTR 120).



02

Integrated control panel

For unfolding the legs.



03

Unfolding

Automatically unfolds the legs while exiting the ambulance.



Comparison of the features of the JUPITER TROLLEY STRETCHER SERIES

| TROLLEY | TG 8802 | TG 8804 |
|---|-------------------|---------|
| 8 positions | ✓ | ✓ |
| Front wheels (200 mm) | Lockable rotating | ✓ |
| | Fixed | ✓ |
| Rear rotating wheels (200 mm) | ✓ | ✓ |
| Trendelenburg and reverse Trendelenburg | ✓ | ✓ |
| Anti-fold safety system trolley legs | ✓ | ✓ |
| Collapsible rear handles | ✓ | ✓ |

| STRETCHER | MTR-120 | MT-118 IN |
|---|---------|-----------|
| Adjustable headrest with shock absorber | ✓ | |
| Fixed leg rest | | |
| Double-jointed leg rest with 2 shock absorbers | ✓ | |
| Longitudinally collapsible guardrails | ✓ | |
| Laterally collapsible guardrails | ✓ | |
| Folding IV pole | ✓ | |
| Aluminium tray for cardiac massage | ✓ | |
| Aluminium tray for incubation equipment | | ✓ |
| 4 telescopic handles | ✓ | ✓ |
| Telescopic central handles | ✓ | |
| Safety belt kit for abdomen and legs | | |
| Safety belt kit for shoulders, abdomen and legs | ✓ | |



Compatibility with rails

| Rails | JUPITER TROLLEY STRETCHER SERIES | JUPITER COMPACT SERIES |
|--------|----------------------------------|------------------------|
| R-419L | ✓ | |
| R-450N | ✓ | |
| R-800L | | ✓ |

Product certification

All KARTSANA products comply with the EN 1865 and 1789 standards.



04

Easy to transfer patients

Laterally and longitudinally collapsible guardrails.



TG-8802

- ▷ 8 intermediary position.
- ▷ Trendelenburg and reverse Trendelenburg.
- ▷ 5 anchoring points to MTR-120 stretchers.



Aluminium



Reclinable



2 rotating wheels



280 kg



TG-8804

- ▷ Same features as TG-8802 except:
4 rotating wheels (200 mm).



Aluminium



Reclinable



4 rotating wheels




280 kg



DIMENSIONS AND WEIGHT

| | TG-8802 | TG-8804 |
|---|-----------|-----------|
| H Height (mm) | 1045 | 1045 |
| HC Height from floor to leading wheel (mm) | 720 | 720 |
| W Width (mm) | 575 | 575 |
| L Minimum/maximum length (mm) | 1990/2369 | 1990/2369 |
| Weight (kg) | 35.9 | 37.9 |



MT-118 IN

- ▷ Versatile, 2-in-1 stretcher.
- ▷ Suitable for all incubator models on the market.



Aluminium



Compatible with incubator



MTR-120

- ▷ Adjustable headrest.
- ▷ Laterally and longitudinally collapsible guardrails.
- ▷ Safety belt kit for shoulders, abdomen and legs.
- ▷ Aluminium tray for performing cardiac massages.
- ▷ Double-jointed leg rest.



Articulated



Guardrails



Aluminium

CERTIFIED
10G



DIMENSIONS AND WEIGHT

| | MT-118 IN | MTR-120 |
|--------------------------------------|-----------|---------|
| H Height (mm) | 150 | 150 |
| W Width (mm) | 590 | 590 |
| L Minimum/maximum length (mm) | 1962 | 1962 |
| Weight (kg) | 32 | 30.6* |



* Weight with mattress and belts.



KARTSANA S.L.

Narcís Monturiol, 34
08192 · Sant Quirze del Vallès
Barcelona (Spain)
T.: +34 93 715 86 72
E.: info@kartsana.com
www.kartsana.com

07-2022

© 2022 · Kartsana. All rights reserved.