

KARTSANA

Life-moving innovation



JUPITER COMPACT SERIES



JUPITER

COMPACT SERIES

JUPITER

Maximum capacity, minimal effort

High-capacity mechanical stretchers that are surprisingly light and manoeuvrable.



01

Load capacity of 280 kg

Stretcher which can lift a load of up to 280 kg unassisted.



02

Easy to transfer patients

Laterally and longitudinally collapsible guardrails.



03

Double-jointed leg rest



04

Bariatric kit



Product certification

All KARTSANA products comply with the EN 1865 and 1789 standards.



COMPACT

Comparison of the features of the JUPITER COMPACT SERIES

		TG-880 C2	TG-880 C4	TG-880S C4	TG-880M C4	TG-880C4 IN
Adjustable headrest with shock absorber		✓	✓	✓	✓	
Double-jointed leg rest with 2 shock absorbers		✓	✓	✓	✓	
Positions	5			✓		
	8	✓	✓		✓	✓
Front wheels (200 mm)	Lockable rotating		✓	✓	✓	✓
	Fixed	✓				
Rear wheels (200 mm)	Lockable rotating	✓	✓	✓	✓	✓
Trendelenburg and reverse Trendelenburg		✓	✓	✓	✓	✓
Longitudinally collapsible guardrails		✓	✓	✓	✓	
Laterally collapsible guardrails		✓	✓	✓	✓	
Folding IV pole		✓	✓	✓	✓	
Safety belt kit for shoulders, abdomen and legs		✓	✓	✓	✓	
Aluminium tray for cardiac massage		✓	✓	✓	✓	
4 telescopic handles		✓	✓	✓	✓	✓

Compatibility with rails

Rails	JUPITER COMPACT SERIES
R-419	Only model TG-880SC4
R-419L	All models except TG-880SC4
R-450N	All models except TG-880SC4
R-450S STAINLESS STEEL	Only model TG-880SC4
R-800G	Only model TG-880SC4
R-800L	All models except TG-880SC4

JUPITER COMPACT Paediatric

The JUPITER SERIES features the TG-880 C4 IN stretcher which is specially designed for paediatric services that require an incubator. More information on page 42.



TG-880 C4 / TG-880 C2

- ▷ 8 intermediary position.
- ▷ Trendelenburg and reverse Trendelenburg.
- ▷ Maximum load 280 kg.



Aluminium



Reclinable



2 rotating wheels



4 lockable wheels



280 kg



CERTIFIED 10G




Paediatric option - TG-880 C4 IN



DIMENSIONS AND WEIGHT

	SIZE L		SIZE M	SIZE S
	TG-880 C2/C4	TG-880 C4 IN	TG-880M C4	TG-880S C4
H Height (mm)	986	986	924	907
HC Height from floor to leading wheel (mm)	720	720	670	601
W Width (mm)	595	595	590	595
L Length (mm)	1974	1974	1990	1975
Weight (kg)	57 / 58.9	55	59,6	59.1







KARTSANA S.L.

Narcís Monturiol, 34
08192 · Sant Quirze del Vallès
Barcelona (Spain)
T.: +34 93 715 86 72
E.: info@kartsana.com
www.kartsana.com

07-2022

© 2022 · Kartsana. All rights reserved.