

# KARTSANA

Life-moving innovation



EVACUATION CHAIRS



# EVACUATION CHAIRS

# CHAIR

## For all needs

Light and safe chairs for transporting and evacuating patients whatever the situation.






## S-127



- ▷ Aluminium folding chair with an epoxy paint finish.
- ▷ 2 safety belts.
- ▷ 2 wheels.




### DIMENSIONS AND WEIGHT

		S-127	S-129	S-129RG
UNFOLDED CHAIR	L Length (mm)	498	720	760
	W Width (mm)	489	500	525
	H Height (mm)	1035	921	990
FOLDED CHAIR	L Length (mm)	498	247	247
	W Width (mm)	196	500	525
	H Height (mm)	596	921	990
Weight (kg)		8.0	8.6	10.9
Load capacity (kg)		150	150	150

## S-129




- ▷ Aluminium folding chair with an epoxy paint finish.
- ▷ 4 wheels.
- ▷ 6 handles (4 collapsible + 2 extendable).
- ▷ Folding footrest.
- ▷ 3 safety belts.
- ▷ Weight: 8,600 kg.



## S-129 RG



- ▷ Aluminium folding chair.
- ▷ 4 wheels.
- ▷ 6 handles (4 collapsible + 2 extendable).
- ▷ 3 safety belts.
- ▷ Folding footrest.
- ▷ Weight: 10,900 kg.





KARTSANA S.L.

Narcís Monturiol, 34  
08192 · Sant Quirze del Vallès  
Barcelona (Spain)  
T.: +34 93 715 86 72  
E.: [info@kartsana.com](mailto:info@kartsana.com)  
[www.kartsana.com](http://www.kartsana.com)

07-2022

© 2022 · Kartsana. All rights reserved.