



DATA SHEET

# S-127 / S-129 / S-129 RG

EVACUATION CHAIRS

MECHANICAL



### S-127:

- ▷ Aluminium folding chair with an epoxy paint finish.
- ▷ 2 safety belts and 2 wheels.

### S-129:


- ▷ Aluminium folding chair with an epoxy paint finish.
- ▷ 4 wheels.
- ▷ 6 handles (4 collapsible + 2 extendable).
- ▷ Folding footrest.
- ▷ 3 safety belts.
- ▷ Weight: 8.600 kg.

### S-129 RG:

- ▷ Aluminium folding chair.
- ▷ 4 wheels.
- ▷ 6 handles (4 collapsible + 2 extendable).
- ▷ 3 safety belts.
- ▷ Folding footrest.
- ▷ Weight: 10.900 kg.

## DIMENSIONS AND WEIGHT

		S-127	S-129	S-129RG
UNFOLDED CHAIR	L Length (mm)	498	720	760
	W Width (mm)	489	500	525
	H Height (mm)	1035	921	990
FOLDED CHAIR	L Length (mm)	498	247	247
	W Width (mm)	196	500	525
	H Height (mm)	596	921	990
Weight (kg)		8.0	8.6	10.9
Load capacity (kg)		150	150	150



**KARTSANA**

Life-moving innovation

Narcís Monturiol, 34  
08192 - Sant Quirze del Vallès  
Barcelona (Spain)

T.: +34 93 715 86 72  
info@kartsana.com  
www.kartsana.com