



DATA SHEET

# TG-880 C4 / TG-880 C2

JUPITER COMPACT SERIES

MECHANICAL



Aluminium



Reclinable



2 rotating wheels



4 lockable wheels





280 kg



- ▷ 8 intermediary position.
- ▷ Trendelenburg and reverse Trendelenburg.
- ▷ Weight 46 kg.
- ▷ Maximum load 280 kg.

## DIMENSIONS AND WEIGHT

	SIZE L		SIZE M	SIZE S
	TG-880 C2/C4	TG-880 C4 IN	TG-880M C4	TG-880S C4
<b>H</b> Height (mm)	986	986	924	907
<b>HC</b> Height from floor to leading wheel (mm)	720	720	670	601
<b>W</b> Width (mm)	595	595	590	595
<b>L</b> Length (mm)	1974	1974	1990	1975
Weight (kg)	57 / 58.9	55	59.6	59.1



## RELATED PRODUCTS



R-450N



R-450S STAINLESS STEEL



TG-880 C4 IN



**KARTSANA**

Life-moving innovation

Narcís Monturiol, 34  
08192 - Sant Quirze del Vallès  
Barcelona (Spain)

T.: +34 93 715 86 72  
info@kartsana.com  
www.kartsana.com