



DATA SHEET

# TG-8802 / TG-8804

## JUPITER TROLLEY STRETCHER

MECHANICAL



Aluminium



Reclinable



2 rotating wheels



4 rotating wheels



280 Kg



CERTIFIED 10G


**\*TG-8804**



- ▷ 8 intermediary position.
- ▷ Trendelenburg and reverse Trendelenburg.
- ▷ 5 anchoring points to MTR-120 stretchers.
- ▷ **TG-8804**: 4 rotating wheels (200 mm).

### DIMENSIONS AND WEIGHT

	TG-8802	TG-8804
<b>H</b> Height (mm)	1045	1045
<b>HC</b> Height from floor to leading wheel (mm)	720	720
<b>W</b> Width (mm)	575	575
<b>L</b> Minimum/maximum length (mm)	1970	1970
Weight (kg)	35.9	37.9



### RELATED PRODUCTS



MT-120



TG-880 C4 IN



**KARTSANA**

Life-moving innovation

Narcís Monturiol, 34  
08192 - Sant Quirze del Vallès  
Barcelona (Spain)

T.: +34 93 715 86 72  
info@kartsana.com  
www.kartsana.com