



DATA SHEET

TG-880 C4 IN

INCUBATOR STRETCHERS

SPECIAL EQUIPMENT



Aluminium



4 lockable wheels



280 Kg



- ▷ Stretcher prepared for pediatric services with incubator.
- ▷ Maximum weight 280 kg.
- ▷ 6 intermediate positions.
- ▷ Trendelenburg and anti-Trendelenburg.
- ▷ 4 telescopic grips.
- ▷ 4 Ø200 mm rotating wheels.




* Upper tray MT-118 IN not included in reference TG-880 C4 IN.

DIMENSIONS AND WEIGHT

	TG-880 C4 IN
H Height (mm)	986
HC Height from floor to leading wheel (mm)	720
W Width (mm)	595
L Length (mm)	1974
Weight (kg)	55

	MT-118 IN*
H Height (mm)	150
W Width (mm)	590
L Length (mm)	1962
Weight (kg)	32



RELATED PRODUCTS



TG-880 C4 / TG-880 C2



MT-118 IN



KARTSANA

Life-moving innovation

Narcís Monturiol, 34
08192 - Sant Quirze del Vallès
Barcelona (Spain)

T.: +34 93 715 86 72
info@kartsana.com
www.kartsana.com