

DATA SHEET
TGE-241

241 COMPACT SERIES

MECHANICAL



Aluminium



Reclinable



250 Kg



Laterally and longitudinally collapsible guardrails



Rotating wheels



CERTIFIED
10G



- ▷ 250kg load capacity
- ▷ Shock-assisted adjustable head
- ▷ Trendelenburg, anti-Trendelenburg and intermediate position
- ▷ 4 Ø200 wheels (2 fixed + 2 rotating with brake)
- ▷ 2 heights available to adapt to the rail ramp of the ambulance (M and XL)
- ▷ Side and longitudinal folding rails
- ▷ Foldable IV pole
- ▷ Shoulder, abdomen and leg belt kit (3 points)
- ▷ Aluminum tray for cardiac massage
- ▷ 4 telescopic grips
- ▷ Automatic leg unfolding
- ▷ Fireproof welded mattress

DIMENSIONS AND WEIGHT

	S	M	L	XL
H Height (mm)	-	860	-	950
HC Height from the floor to the leading wheel (mm)	-	585	-	700
W Width (mm)	-	590	-	590
L Length (mm)	-	1951	-	1951
Weight (kg)	-	46.5	-	44.1



RELATED PRODUCTS



R-419



KARTSANA

Life-moving innovation

Narcís Monturiol, 34
08192 · Sant Quirze del Vallès
Barcelona (Spain)

T.: +34 93 715 86 72
info@kartsana.com
www.kartsana.com