

DATA SHEET
TGX-241

241 COMPACT SERIES

MECHANICAL



Aluminium



Reclinable



250 Kg



Laterally and longitudinally collapsible guardrails



Leg rest with shock absorber



Foldable


CERTIFIED
10G



- ▷ 250kg load capacity
- ▷ Shock-assisted adjustable head
- ▷ Trendelenburg, anti-Trendelenburg and intermediate position
- ▷ 4 Ø200 wheels (2 fixed + 2 rotating with brake)
- ▷ 3 heights available to adapt to the ramp of the rail of the ambulance (M, L and XL)
- ▷ Side and longitudinal folding rails
- ▷ Foldable IV pole
- ▷ Shoulder, abdomen and leg belt kit (3 points)
- ▷ Aluminum tray for cardiac massage
- ▷ 4 telescopic grips
- ▷ Automatic leg unfolding
- ▷ Trigger anti-accidental actuation safety
- ▷ feet (TGX-241M only)
- ▷ Folding head with collapsible grips for maneuvering in hard-to-reach places
- ▷ Shock-assisted single-joint leg
- ▷ Fireproof welded mattress

DIMENSIONS AND WEIGHT

| | S | M | L | XL |
|---|---|------|------|------|
| H Height (mm) | - | 860 | 909 | 953 |
| HC Height from floor to leading wheel (mm) | - | 590 | 642 | 686 |
| W Width (mm) | - | 590 | 590 | 590 |
| L Length (mm) | - | 1953 | 1953 | 1953 |
| L2 Length folded (mm) | - | 1630 | 1630 | 1571 |
| Weight (kg) | - | 49.5 | 51.5 | 53 |



RELATED PRODUCTS



R-419



KARTSANA

Life-moving innovation

Narcís Monturiol, 34
08192 - Sant Quirze del Vallès
Barcelona (Spain)

T.: +34 93 715 86 72
info@kartsana.com
www.kartsana.com