



DATA SHEET

TG-1200 SILVER+

SILVER+ SERIES

ELECTRIC EQUIPMENT



Aluminium



Electric



250 Kg




CERTIFIED 10G



- ▷ Control of the stretcher from the technician's position.
- ▷ Load capacity of 250 kg.
- ▷ Manual motor release in case of emergency.
- ▷ Individual battery with extra capacity of 50 cycles, tested with 100 Kg of weight.
- ▷ Adjustable head assisted by shock absorber.
- ▷ Shock-absorber-assisted single-joint leg.
- ▷ 4 rotating wheels (2 with brake + 2 lockable).
- ▷ Longitudinal and lateral folding handrails.
- ▷ Folding IV holder.
- ▷ Three-point belt kit.
- ▷ Aluminum tray for cardiac massage.

PRODUCT SPECIFICATIONS

	TG-1200
H Max. height (mm)	1135
HC Max. adjustable height (mm)	786
H2 Min. height (mm)	350
L Length (mm)	1955
W Width (mm)	566
Load capacity (Kg)	230
Weight (Kg)	73



RELATED PRODUCTS



R-1200 SILVER+



INSTRUMENTAL TABLE C-101



KARTSANA

Life-moving innovation

Narcís Monturiol, 34
08192 · Sant Quirze del Vallès
Barcelona (Spain)

T.: +34 93 715 86 72
info@kartsana.com
www.kartsana.com

04-2024